drive

DeVilbiss®
HEALTHCARE

1060AW

en Drive DeVilbiss® PulmO2™ 10-Liter Oxygen Concentrator Instruction Guide

WARNING—Read instruction guide before operating this equipment.

ASSEMBLED IN USA

 **NO SMOKING**

Rx
ONLY

Scan here for Video
Instructions

es Drive DeVilbiss® PulmO2™ 10-Liter Oxygen Concentrator Instruction Guide

WARNING—Read instruction guide before operating this equipment.

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 **INTERDICTION DE FUMER**

Rx
ONLY

NOTE: NOT MANUFACTURED WITH NATURAL RUBBER LATEX

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Read this entire guide before using your PULMO2™ oxygen concentrator. Important safeguards are indicated throughout this guide.

WARNING

Under certain circumstances, oxygen therapy can be hazardous. Seek medical advice before using an oxygen concentrator.

1. WARNINGS AND PRECAUTIONS

Pay special attention to all safety information. Imminently and potentially hazardous information is highlighted by these terms:

DANGER

Indicates an imminently hazardous situation which could result in death or serious injury to the user or operator if not avoided.

WARNING

Indicates a potentially hazardous situation which could result in death or serious injury to the user or operator if not avoided.

CAUTION

Indicates a potentially hazardous situation which could result in property damage, injury, or device damage if not avoided.

IMPORTANT

Indicates important information you should know.

NOTE

Indicates notes, useful tips, recommendations, and information.

READ ALL INSTRUCTIONS BEFORE USE.

DANGER

NO SMOKING signs should be prominently displayed.

- ▶ Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy.
- ▶ Do not leave the nasal cannula or mask on bed coverings or chair cushions if the oxygen concentrator is turned on but not in use. The oxygen will make the materials flammable. Turn the oxygen concentrator OFF when not in use.
- ▶ Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located.
- ▶ If you intend to smoke, you must always turn the oxygen concentrator OFF, remove the cannula and leave the room where either the cannula or mask or the oxygen concentrator is located. If unable to leave the room, you must wait 10 minutes after you have turned off the oxygen concentrator before smoking.
- ▶ Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.
- ▶ Open flames during oxygen therapy are dangerous and are likely to result in fire or death. Do not allow open flames near any oxygen carrying accessories.
- ▶ Drive DeVilbiss oxygen concentrators are equipped with a fire mitigating outlet fitting that prevents propagation of fire into the unit.

WARNINGS

- To avoid electric shock, do not plug the concentrator into an AC outlet if the concentrator cabinet is broken. Do not remove the concentrator cabinet. The cabinet should only be removed by an appropriately trained service technician. Do not apply liquid directly to the cabinet or utilize any petroleum-based solvents or cleaning agents.
- Improper use of the power cord and plugs can cause a burn, fire or other electric shock hazards. Do not use the unit if the power cord is damaged.
- Ensure the power cord is fully inserted into the concentrator connector and a fully functioning AC wall outlet.
- Failure to do so may cause an electrical safety hazard.
- The accessories (nasal cannula, masks, oxygen tubing, humidifiers, etc.) that supply oxygen to the patient must be equipped with a means that, in case of fire, stops the propagation of fire through the accessory for the safety of the patient and others. A fire activated flow-stop or thermal fuse device, if available, should be used with oxygen supply accessories. These types of flow-stop devices stop the flow of oxygen to the patient in the event of fire. This means of fire protection should be located as close to the patient as practicable.
- Locate oxygen tubing and power supply cords to prevent tripping hazards and reduce the possibility of entanglement or strangulation.
- Do not lubricate fittings, connections, tubing or other accessories of the oxygen concentrator to avoid the risk of fire and burns.
- Do NOT use lubricants, oils or grease during use or service due to the risk of fire.
- Before attempting any cleaning procedures, turn the unit OFF.
- Use only water-based lotions or salves that are oxygen-compatible before and during oxygen therapy. Never use petroleum-based or oil-based lotions or salves to avoid the risk of fire and burns.
- Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.
- If you feel discomfort or are experiencing a medical emergency while undergoing oxygen therapy, seek medical assistance immediately to avoid harm.
- Geriatric, pediatric or any other patient unable to communicate discomfort can require additional monitoring and/or a distributed alarm system to convey the information about the discomfort and/or the medical urgency to the responsible caregiver to avoid harm.
- Use of this device at an altitude above 6,562 feet (2000 meters) or above a temperature of 104°F (40°C) or greater than 95% relative humidity may affect the flow rate and the percentage of oxygen and consequently the quality of the therapy. Refer to specifications for details regarding parameters tested.
- To ensure you receive the therapeutic amount of oxygen delivery according to your medical condition, the Oxygen Concentrator must:
 - a. be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels.
 - b. be used with the specific combination of parts and accessories that are in line with the specification of the concentrator manufacturer and that were used while your settings were determined.
- Your delivery settings of the oxygen concentrator should be periodically reassessed for the effectiveness of therapy.
- For your safety, the oxygen concentrator must be used according to the prescription determined by your physician.
- Under certain circumstances, oxygen therapy can be hazardous. Seek medical advice before using an oxygen concentrator.
- The proper placement and positioning of the prongs of the nasal cannula in the nose is critical to the amount of oxygen delivered to the respiratory system of the patient.
- Do not modify the oxygen concentrator or accessories.
- Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- This device may be affected by children unintentionally changing the device setting or inserting foreign objects into device openings. Supervise children when the device is in operation.
- The device can be affected by pet hair or insects getting into the device or obstructing cooling vents on the device.
- Concentrated oxygen flow rates above 1 LPM for Neonatal and 2 LPM for infant patients can create a risk of high FiO₂ (>50%) which may cause issues like retinopathy of prematurity and bronchopulmonary dysplasia. In infants, very high FiO₂s may lead to oxygen toxicity - damage to lung tissue.

CAUTION

Risk of injury or damage

- Operate the unit in a cool, dry area with good ventilation, located on a hard surface, avoid thick rugs or carpeting. NEVER block the air intake or exhaust vents
- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Use only recommended cleaning and disinfection agents to avoid damage to the device.
- Keep exposed body parts, such as hands and feet, a minimum of 46 inches (1.2 meters) away from the exhaust vents to avoid the risk of burns. Single fault conditions may result in visual and audible alerts and alarms. The oxygen concentrator should be located in a well-ventilated area. DO NOT operate the unit in a closed or confined space, such as a closet, bathroom, etc. Avoid operating the device near smoke pollutants and fumes.
- To prevent product damage, do not attempt to operate the unit without the cabinet air filter or while the filter is still damp.
- Incompatible parts or accessories can result in degraded performance. Use only recommended parts and filters to ensure reliable operation of the product.
- To prevent product damage, do not attempt to operate the unit without the compressor intake filter.
- Use of harsh chemicals (including alcohol) is not recommended. If bactericidal cleaning is required, a non-alcohol based product should be used to avoid inadvertent damage
- Install, use and maintain the PulmO2 10 Liter Oxygen Concentrator following this instruction guide to minimize the environmental impact of the oxygen concentrator through its expected life. The oxygen concentrator consumes maximum 520 watts during normal use. Water can be used as needed for humidification of the delivered oxygen. The Compressor Intake Filter and Final Hepa Filter are disposable and could need to be replaced twice during the expected service life. Other byproducts created during the normal use of the oxygen concentrator are heat, noise and nitrogen gas. Follow instructions to minimize the effects of heat and noise. Nitrogen gas discharged by the oxygen concentrator is quickly dissipated back into room air. The oxygen concentrator contains no hazardous substances.

NOTE

- Under normal and single fault conditions, the concentrator releases warm air out the bottom of the unit (exhaust vents) which may discolor temperature sensitive flooring surfaces. The concentrator should not be used over flooring that is sensitive to heat staining. The Manufacturer is not responsible for flooring that becomes discolored.
- Under extreme environmental conditions and a single fault condition occurs, the following device surface temperatures may exceed 106 °F (41 °C). See Table 1 below for the model specific maximum temperature and safe contact guidance:

CAUTION

Table 1 - 1060AW

Under extreme environmental conditions and a single fault condition occurs, the following device surface temperatures may exceed 106 °F (41 °C). See Table 1 below for the model specific maximum temperature and safe contact guidance:

Description	Maximum Temperature		Max safe contact time
	1060AW		
	°F	°C	
Outlet Vents	145.8	63.2	Less than 1 minute
Oxygen outlet fitting	107.6	42.0	Less than 10 minutes
Oxygen Tubing Coupler	106.0	41.1	Less than 10 minutes
Cannula tubing at the outlet fitting	106.2	41.2	Less than 10 minutes

CAUTION

To avoid possible injury or loss of therapy, turn the device off and disconnect from patient before moving the device.

SAVE THESE INSTRUCTIONS.

1.1 INTRODUCTION

This instruction guide will acquaint you with device. Make sure that you read and understand this guide before operating your unit. Important safeguards are indicated throughout this guide. Pay special attention to all safety information. Contact your equipment provider should you have any questions.

1.2 Intended Use

The PulmO2 10 Liter oxygen concentrator is intended to provide supplemental oxygen to patients requiring oxygen therapy.

1.3 Contraindications

None known.

1.4 Indications For Use

The PulmO2 10 Liter oxygen concentrator is intended to provide supplemental oxygen to patients requiring oxygen therapy. The device may be used in the home or a clinical setting. The device is not intended to be life sustaining or life supporting.

1.5 Essential Performance

The essential performance of this device is to deliver a continuous flow of oxygen. Visual and audible alarms are provided if the flow of the of the concentrator is lost or degraded.

1.6 About Your Device

Oxygen concentrators are the most reliable, efficient and convenient sources of supplemental oxygen available today. The oxygen concentrator is electrically operated. The unit separates oxygen from room air which allows high-purity supplemental oxygen to be delivered to you through the oxygen outlet. Although the concentrator filters the oxygen in a room, it will not affect the normal amount of oxygen in your room.

NOTE

- ▶ Drive DeVilbiss recommends for optimal service life that the PulmO2 10 Liter Oxygen Concentrator be operated for at least 30 minutes after it is powered on. Shorter periods of operation, operating in extreme temperature humidity conditions or in the presence of contaminants, and/or handling and storage conditions outside those specified, may affect the long term reliable operation of the product.

DRIVE DEVILBISS (OSD)[®] OXYGEN SENSING DEVICE

The OSD (Oxygen Sensing Device) is a device within your concentrator that monitors the oxygen produced by your unit.

The OSD lights on the top panel are defined as follows:

- ▶ Green Good Oxygen LED – acceptable oxygen level.
- ▶ Yellow Low Oxygen LED – below the good oxygen level.
- ▶ Red Critical Low Oxygen LED - below the critical oxygen level.

If the oxygen purity falls below the acceptable level: The green Good Oxygen LED will shut off, the yellow Low Oxygen LED will light up, and an intermittent audible signal will sound. If the oxygen purity continues to fall into the critical low level, the Red Critical Low Oxygen LED alerts with an audible alarm.

Refer to the Troubleshooting (**Section 4**) in this guide. Do not attempt any other maintenance. Contact your service technician immediately if the issue is not resolved.

IMPORTANT

- ▶ It is recommended that the clinical staff or provider lock the flow control knob to prevent inadvertent adjustment. A flow setting other than prescribed may affect the patient therapy.
- ▶ Do not service or clean this device while in use with a User.
- ▶ The 1060AW Series oxygen concentrator, its parts and accessories are specified for use at flow rates between 1 LPM and 10 LPM.
- ▶ When used in a homecare environment, the patient is the intended operator.

1.7 Important Parts of Your Concentrator


Scan for video Instructions

Please take time to familiarize yourself with your PulmO2 10 Liter oxygen concentrator before operating.

Front & Back View (Figure A & B)

1. Control Panel / Operation (LED indicator panel)
2. Power Switch
| = ON
O = OFF
3. Flow meter knob
4. Flow meter
5. Oxygen outlet – oxygen is dispersed through this port
6. Handle (top) and Handgrips (both sides)
7. Filter Door
8. Compressor Intake Filter
9. Hour meter
10. Exhaust Vents
11. IEC power connector
12. Line cord strap (not shown)
13. Cabinet Air Filter (gross particle)
14. (2 locking, 2 non-locking) Casters
15. Air Intake Vent
16. Filter Door Vents (behind cabinet air filter #13)
17. Oxygen Outlet Tubing Connector CN 100





⚠ WARNING

When the device is used under extreme operating or single fault conditions, the exhaust air near the exhaust vents on the bottom of the unit may exceed 41° C.

Control Panel Label (Figure C & D)

NOTE - Symbols on control panel are outlined b&w when not active and illuminated/backlit when active.

1. O₂ % (Oxygen Percentage) Status LEDs (see page EN-16)
 - a. Startup Low O₂ yellow LED
 - b. Critical Low Oxygen red LED
 - c. Low Oxygen yellow LED
 - d. Good Oxygen green LED
2. Clinical Staff Recurring Reminders (see page EN-16)
 - a. Check Cabinet Air Filter LED - Illuminates weekly (170 hours), auto-resets after 12 hours.
 - b. Check Compressor Intake Filter LED - Illuminates every 3 months (3200 hours), auto-resets after 12 hours.
3. Clinical Staff and Service Technician Attention Required Alarms
 - a. Low Flow Condition - check for condition cause when illuminated
 - b. Over Temperature Condition - check for condition cause when illuminated
 - c. Service Required - if LED activates, discontinue use and notify service technician to repair unit.
4. Alarm Acknowledge and Active LED - Silences audible alarm and illuminates when pressed/active.

1.8 ACCESSORY EQUIPMENT

Use of Accessories

CAUTION

Use only recommended parts and filters to ensure reliable operation of the product.

Recommended Accessories

High Flow (6-15 LPM) Bubble Humidifier..... Salter Labs 7900 or equivalent
Low Flow (up to 6 LPM) Bubble Humidifier Salter Labs 7600 or equivalent
High Flow Nasal Cannula..... Salter Labs 1600HF or equivalent

Refer to the accessory manufacturer's instructions for use. A mask or any nasal cannula can be used with continuous flow delivery and may be sized according to your prescription as recommended by the clinical staff who should also give you advice on the proper usage, maintenance, and cleaning.

WARNING








- ▶ The accessories (nasal cannula, masks, oxygen tubing, humidifiers, etc.) that supply oxygen to the patient must be equipped with a means that, in case of fire, stops the propagation of fire through the accessory for the safety of the patient and others. A fire activated flow-stop or thermal fuse device, if available, should be used with the oxygen supply accessories. These types of flow-stop devices stop the flow of oxygen to the patient in the event of fire. This means of fire protection should be located as close to the patient as practicable.
- ▶ Do not use accessories, detachable parts and materials not described in the instructions for use. Do not modify the oxygen concentrator or accessories.
- ▶ It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen – consult your physician.
- ▶ There are many types of humidifiers, oxygen tubing and cannulas/masks that can be used with this device. Certain humidifiers and accessories may impair the device's performance.
- ▶ Do not interconnect this oxygen concentrator with other equipment not described in the instructions for use.



NOTE

- The device should only be used with bubble humidifiers that are designed for use with flows up to 10 liters per minute and 8.5 psi pressure.
- A maximum of 50 feet (15 meters) of crush-proof oxygen tubing, plus 7 feet (2.1 meters) of cannula, plus a bubble humidifier is allowed between the concentrator and the patient.
- The clinical staff or your provider should verify the compatibility of the oxygen concentrator and all of the parts used to connect to the patient before use.
- The clinical staff or your provider may have preset the flow meter so that it cannot be adjusted.
- Your physician has prescribed either a nasal cannula or face mask. In most cases, they are already attached to the oxygen tubing. If not, follow the manufacturer's instructions for attachment.

1.9 SYMBOLS GLOSSARY



	Manufacturer		Medical Device		Serial Number		LOT Number		Clinical Staff
	Catalog Number		Class II Equipment		Attention Required		Recurring Reminders		Service Technician
	Humidifier Bottle		Place on a hard surface		Well ventilated area		30 cm away from wall		On Off
	Operating Temperature Range +5 to +40°C (+41 to +104°F)		Operating Atmospheric Pressure Range 795 to 1013 hPa (Sea level to 2000m (6562 ft))		TUV Rheinland C-US approval mark		Type BF applied part		User
	Storage Humidity Limit		Keep Dry		Storage Temperature Limits		Hour Meter		Alternating Current
	MRI Unsafe symbol – Unsafe for Magnetic Resonance Environment <i>This symbol has a red circle and diagonal bar on the product label.</i>		Do not use wet cabinet air filter <i>This symbol has a red circle and diagonal bar on the product label.</i>		Do not cover <i>This symbol has a red circle and diagonal bar on the product label.</i>		Do not Disassemble <i>This symbol has a red circle and diagonal bar on the product label.</i>	 DO NOT USE NEAR SURGICAL EQUIPMENT with High Frequency Surgical Equipment / Radio Frequency Equipment <i>This symbol has a red circle and diagonal bar on the product label.</i>	
	CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.		Over Temperature LED <i>This symbol illuminates Red when active</i>		Startup Low O2 Indicator <i>Oxygen is low during initial start up. This symbol illuminates Yellow when active.</i>		Alarm Acknowledge Button and Active LED <i>This symbol illuminates White when active.</i>	 DO NOT USE NEAR SURGICAL EQUIPMENT with High Frequency Surgical Equipment / Radio Frequency Equipment <i>This symbol has a red circle and diagonal bar on the product label.</i>	
	Check Cabinet Air Filter Notification LED <i>This symbol illuminates Yellow when active</i>		Check Compressor Intake Filter Notification LED <i>This symbol illuminates Yellow when active</i>		Low Flow LED <i>This symbol illuminates Red when active</i>		Service Required LED <i>This symbol illuminates Red when active.</i>		Oxygen Outlet marked with rated range of gas pressure and rated range of gas flowrate
	It is mandatory to read and understand the operating instructions prior to use. <i>This symbol has a blue background on the product label.</i>		Electric Shock Hazard. Cabinet to be removed by authorized personnel only. <i>This symbol has a yellow background on the product label.</i>		No smoking near patient or device. <i>This symbol has a red circle and diagonal bar on the product label.</i>		Use no Oil, Grease or Lubricants. <i>This symbol has a red circle and diagonal bar on the product label.</i>		Do not use near heat or open flames <i>This symbol has a red circle and diagonal bar on the product label.</i>

	<p>This device contain electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU- Waste Electrical and Electronic Equipment (WEEE)</p>		<p>General Warning <i>This symbol is used throughout this manual to indicate hazardous situations to avoid.</i></p>		<p>Important Information <i>This symbol is used throughout this manual to indicate important information you should know.</i></p>		<p>Note and Information Symbol <i>This symbol is used throughout this manual to indicate notes, useful tips, recommendations and information.</i></p>		<p>Ingress Protection - Protected against finger access to hazardous parts; protected against vertically falling water drops.</p>		
 <table border="1" data-bbox="162 430 470 546"> <tr> <td>Critical Low Oxygen ● Illuminates Red when active</td> <td>Low Oxygen LED ● Illuminates Yellow when active</td> <td>Good Oxygen LED ● Illuminates Green when active</td> </tr> </table> <p>O2 % (Oxygen Percentage) Status LEDs</p>			Critical Low Oxygen ● Illuminates Red when active	Low Oxygen LED ● Illuminates Yellow when active	Good Oxygen LED ● Illuminates Green when active		<p>Maximum recommended flow rate: 10 LPM</p>		<p>QR code link to user and service information</p>		
Critical Low Oxygen ● Illuminates Red when active	Low Oxygen LED ● Illuminates Yellow when active	Good Oxygen LED ● Illuminates Green when active									



HOW TO SET UP THE DEVICE

1. With the power switch (Figure A) of the device in the OFF position, insert the power cord (Figure B) into the IEC Connector located on the backside of the device. Insert the other end of the power cord into the wall outlet as shown below. (Figure C)

2. With a clean cloth wipe away any dust that may be on the cabinet air filter (Figure D) and filter door vents (Figure E), which are located on the back of the device, as shown below.

3. Connect oxygen tubing to the oxygen outlet by either installing a humidifier bottle or attaching the tubing directly.


3.1 Humidifier:

- a. Fill the humidifier bottle as per manufacturer's instructions.
- b. Tightly screw the wing nut located on the top of the humidifier bottle to the oxygen outlet so that the humidifier bottle hangs from the oxygen outlet, as shown. (Figure F)
- c. Attach the oxygen tubing directly to the humidifier bottle outlet fitting, as shown. (Figure G)

3.2 No Humidifier:

- a. Attach Oxygen Outlet Connector to the oxygen outlet
- b. Attach the Oxygen Tubing to the Outlet Connector. (Figure H)


HOW TO USE THE DEVICE

1. With the oxygen tubing connected to either the humidifier OR the oxygen outlet connector, press the power switch of the unit to the **“ON”** position.

When the Oxygen Concentrator is first turned **“ON”**, you will see all ten (10) LEDs on the control panel briefly light up and you should hear an alarm. The Oxygen Concentrator will then operate in **“start-up”** mode for up to **15 minutes** with the Startup Low O₂ LED light on until a normal oxygen level is achieved, at which time the Green Good O₂ LED lights up and will remain lit.

2. Check that the flow meter is set to your prescribed oxygen flow rate. Your oxygen concentrator is now ready for use.
3. If using an oxygen cannula (see Figure I), insert the nasal prongs of your oxygen cannula into nose, OR if using an oxygen mask (see Figure J) position the oxygen mask over your nose and mouth. Wrap the tubing of your oxygen cannula/mask over your ears and place in front of your body. With the oxygen cannula/mask secured, breathe normally to receive oxygen therapy.


4. When finished using the unit, press Power Switch to turn unit **OFF**.

The healthcare provider is to assess the needs of the patient, at time of installation, for backup supplies of supplemental oxygen in case of oxygen concentrator or power failure, based on the condition of the patient, the environment in which the patient lives and the ability to resupply the patient with backup supplies of supplementary oxygen. The healthcare provider should reassess as conditions change.

If the power cord needs to be removed.

1. Press the yellow release tab
2. Pull the plug while holding the release tab

WARNINGS

As part of the normal start-up process, all ten (10) LEDs on the control panel should light up and an alarm should sound when the unit is first turned ON. If ANY front panel LEDs DO NOT light up or the alarm DOES NOT sound, the alert system is not functioning properly. Refer to the Troubleshooting chart and contact your service technician if necessary.

- ▶ Ensure the power cord is fully inserted into the oxygen concentrator connector and the power cord plug is completely inserted into a fully functioning AC wall outlet. Failure to do so may cause an electrical safety hazard.
- ▶ Improper use of the power cord and plugs can cause a burn, fire or other electric shock hazards. Do not use the unit if the power cord is damaged

CAUTION

It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen – consult your physician.

NOTE

- ▶ Do not connect to an electrical outlet controlled by a wall switch. If using a switched outlet, make sure the outlet is always switched on.
- ▶ Your healthcare provider or clinical staff should verify the compatibility of the oxygen concentrator and all of the parts used to connect to the patient before use.
- ▶ For prescriptions of 10 LPM, be sure the ball is centered on the 10 liter line. The ball should not touch the red line. Setting the flow higher than 10 may cause the oxygen purity level to drop.
- ▶ The low-flow alarm may activate if the flow meter ball is set below 1 lpm. The unit will continue to run; however, the Service Required light will come on accompanied by an audible alarm. Adjust the flow meter to your prescribed flow.
- ▶ To check your oxygen concentrator and accessories for proper operation; 1. Check the output flow by placing the end of the nasal cannula under the surface of a half-full cup of water and look for the bubbles. 2. Check the system for leaks by bending the nasal prongs over and squeeze tight to stop the flow of oxygen. Look at the flow meter to see that the indicator ball on the flow meter drops to zero. If the indicator ball does not drop to zero, check all connections for possible leaks. Parts to check for leaks are: tubing connections, humidifier bottle and other accessories like firebreaks. Repeat these steps until the flow meter ball drops to zero. Contact your provider or service supplier immediately if you encounter any problems.
- ▶ If the audible signal alarms but the unit is not putting out oxygen, there is no power to the unit. Refer to the Troubleshooting chart on page 10, and contact your service technician if necessary.
- ▶ If an audible low-frequency vibration sound is detected, the unit is not operating properly. Refer to the Troubleshooting chart on page 10, and contact your service technician if necessary.
- ▶ Operate the unit in a cool, dry area with good ventilation, located on a hard surface. Avoid thick rugs or carpeting. NEVER block the air intake or exhaust vents. Keep the unit a minimum of 12 inches (30.5 cm) away from any wall, draperies, or any other objects that might prevent the proper flow of air in and out of your oxygen concentrator. Proper air flow is needed to prevent overheating of the oxygen concentrator. DO NOT place the concentrator near any heat source such as hot air registers or heaters. Overheating of the oxygen concentrator may lead to low oxygen output and a risk of burns. DO NOT operate the unit in a closed or confined space, such as a closet, bathroom, etc. Avoid operating the unit near smoke pollutants and fumes.
- ▶ Your physician has prescribed either a nasal cannula or face mask. In most cases, they are already attached to the oxygen tubing. If not, follow the manufacturer's instructions for attachment.
- ▶ Your clinical staff should assess the needs of the patient for backup supplies for supplementary oxygen in the event of device or power failure.

3. CLEANING AND MAINTENANCE

Scan for video Instructions

ROUTINE CLEANING

• CANNULA/MASK, TUBING AND HUMIDIFIER BOTTLE

Clean and replace the cannula/mask, tubing and humidifier bottle according to the manufacturer's instructions.

• OXYGEN OUTLET CONNECTOR

Remove the oxygen outlet connector and clean as needed or replace between patients.

• FILTER DOOR VENTS

Inspect the vents periodically, and wipe with a dry cloth as needed to remove dust

• CABINET AIR FILTER

The cabinet air filter should be inspected when the Check Cabinet Filter LED illuminates and cleaned as needed by the user or caregiver. Replace if torn or damaged. To clean, follow these steps:

1. Remove the cabinet air filter located on the back of the unit.
2. Wash in a solution of warm water and dishwashing detergent. Rinse thoroughly with warm tap water and towel dry. The filter should be completely dry before reinstalling.

• COMPRESSOR INTAKE FILTER

The compressor intake filter should be inspected when the Check Compressor Intake Filter LED illuminates and replaced as necessary. To check or replace, follow these steps:

1. Open the filter door on the rear of the concentrator.
2. Remove the filter from behind the door and inspect the color and debris. If filter is discolored to a dark gray, it should be replaced.

• EXHAUST VENTS (BOTH SIDES)

Inspect exhaust vents and clean to remove dust or debris with a small brush or cloth as needed. The exhaust vents must be clear of dirt and debris to prevent overheating of the device.

EXTERIOR CABINET

Clean the concentrator exterior cabinet weekly by using a damp cloth and wiping dry; the vents can also be wiped with a damp cloth.

Cleaning

	Recommended cleaning interval	Number of cleaning cycles *	Compatible cleaning method
Outer Cabinet	7 days	260	Water, use only a damp cloth
Filter Door Vents	7 days	260	Wipe with dry cloth, or a cloth dampened with water to remove dust.
Oxygen Outlet Connector	7 days	104	Mild dish soap (2 tbsp)(30ml) and warm water (2 cups)(400 ml)
Cabinet Air Filter	7 days inspect, clean or replace as needed	104	Mild dish soap (2 tbsp) (30ml) and warm water (2 cups) (400 ml), towel dry
Compressor Intake Filter	14 days inspect, replace as needed	---	Do not clean, replace when discolored to dark gray.

* number of cleaning cycles over the expected service life

DANGER

- Before attempting any cleaning procedures, turn the unit OFF. Risk of electrical shock.

3.1 CLEANING AND DISINFECTION FOR SERVICE TECHNICIANS

1. Use only recommended cleaning and disinfection agents.
2. The disinfection process can only be completed by the manufacturer or by an appropriately trained individual. Always wear personal protective equipment when performing this procedure. Use suitable gloves and safety glasses. Cover exposed skin on arms to prevent accidental contact with Sodium Hypochlorite solution that has been applied to the concentrator.
3. Dispose of all accessories that are not suitable for reuse. This includes but may not be limited to the oxygen tubing, tubing connectors, nasal cannula and/or mask, oxygen outlet connector, and humidifier bottle. Reusing accessories marked for single use may cause risk of cross contamination or loss of performance.
4. Clean the exterior of the concentrator with a clean lint-free cloth. Heavy soil should be removed with a clean lint-free cloth dampened with water. A soft bristled brush dampened with water can be used to remove stubborn soil. Dry the concentrator using a clean lint-free cloth if water was used to remove soil.
5. Use 5.25% chlorine Sodium Hypochlorite (Clorox Regular Liquid Bleach or equivalent). Mix one (1) part Sodium Hypochlorite with four (4) parts water in an appropriate clean container. This ratio produces a one (1) part Sodium Hypochlorite to five (5) total parts solution (1:5). The total volume (amount) of solution required is determined by the number of concentrators in need of disinfection.
6. Apply the Sodium Hypochlorite solution in an even manner to the cabinet and power cord using a clean lint-free cloth. The cloth should be dampened only and not dripping of solution. Do not use a spray bottle to apply the solution. Do not saturate the device with the solution. Take care that no solution enters the vent areas on the concentrator base or the area on the back of the unit. Avoid over-saturating the cabinet seams so that no solution residue builds up in these areas. Avoid the caster wells located on the bottom of the unit.
7. Exposure time of the disinfectant solution should be 10 minutes minimum to 15 minutes maximum.
8. After the recommended exposure time, all surfaces of the concentrator should be wiped with a clean lint-free cloth dampened with drinking quality water no warmer than room temperature. Dry the unit with a dry, clean lint-free cloth. This is to remove residue that may stain or leave a film on the unit, especially after repeated disinfections.
9. Check the cord, the plug on the back of the device, the power switch, and the indicator LEDs for possible damage. Replace all damaged or worn components.
10. Inspect Cabinet Air Filter and wash as needed. Replace if torn or damaged.
11. Check the oxygen concentration. If the device is within specification, the compressor intake filter does not need to be replaced between patients. If the oxygen concentration is not within specification, the provider should refer to the service manual section on Troubleshooting.

NOTE

- ▶ There is no portion of the gas pathways through the concentrator that could be contaminated with body fluids under normal conditions. The device patient connection may unintentionally become contaminated with expired gases for a single fault condition i.e., a hose internal to the device becomes disconnected. This condition will cause no flow out of the device and/or an alarm condition. Should this occur, refer to the service manual for additional instructions.
- ▶ Drive DeVilbiss Healthcare recommends that at least the following procedures be carried out by the manufacturer or a qualified third party between uses by different patients.
- ▶ The disinfection process can only be completed by the manufacturer or by an appropriately trained individual.
- ▶ This is a suggested maintenance and service schedule for home oxygen clinical staffs. Individual maintenance requirements may vary based upon local operating conditions, regulations, or other circumstances

DISINFECTION

	Recommended disinfection interval	Number of disinfection cycles	Compatible disinfection method
Cabinet, power cord	Once per week when used with multiple patients.	260	1:5 chlorine Sodium Hypochlorite (5.25%)
Oxygen tubing, tubing connectors, nasal cannula/mask, oxygen outlet connector, humidifier bottle, cabinet air filter	Do not disinfect, replace between patients	N/A	N/A

NOTE

Frequency of inspection and cleaning of filter may be dependent upon environmental conditions like dust and lint.

CAUTION

Do not exceed the recommended number of cleaning/disinfection cycles, device damage may occur.

3.2 SERVICE AND MAINTENANCE FOR SERVICE TECHNICIANS

Recommendations for preventative maintenance at 3-year intervals are outlined in the Service and Maintenance Guidance section below.

PULMO2 10 Liter Oxygen Concentrator Preventive Maintenance/Service Guide

Model	Oxygen Purity Verification	Compressor Intake Filter	Internal Compressor Filter **	Final HEPA Filter **	Cabinet Air Filter	Sieve Beds **	Compressor Cup Seals **
1060AW	Every 3-years or between patients, whichever comes first or with multiple patients, every 6 months	Inspect when Recurring Reminder illuminates on control panel. Replace if filter is dirty or damaged.	Inspect in conjunction with compressor service. Replace if filter is dirty or damaged.	Inspect in conjunction with compressor service. Replace if filter is dirty or damaged.	Inspect when Recurring Reminder illuminates on control panel and wash as needed. Replace if torn or damaged.	When indicated by device performance below specification for oxygen purity, operating pressures and/or other indications of component wear	When indicated by device performance below specification for oxygen purity, operating pressures and/or other indications of component wear

** Sieve bed, compressor cup seals, internal compressor filter and final HEPA filter service should only be performed by appropriately trained service technician.

NOTE

- This is a suggested maintenance and service schedule for home oxygen providers. Individual maintenance requirements may vary based upon local operating conditions, regulations, or other circumstances
- The Manufacturer will make available a Service Manual LT-2355 on request which contains information to assist service personnel with repairing files servicable components.



4. ALARMS AND TROUBLESHOOTING

The following troubleshooting chart will help you analyze and correct minor oxygen concentrator malfunctions. If the suggested procedures do not help, call your service technician. Do not attempt any other maintenance.

This device contains an alarm system which monitors the state of the device and alerts of abnormal operation, loss of essential performance or failures. Alarm conditions are shown on the control panel. The alarm system functions are tested at power up by lighting all visual alarm indicators and sounding the audible alarm (beep).

All alarms are Low Priority Technical Alarms.

Alert or Alarm Condition	LED Icon	Details of Alert or Alarm Condition	Visual Alert or Alarm	Audible Alarm	Action
Power-On Test	ALL	Unit was just turned On; test of all visual/audible indicators	All LEDs ON for a few seconds	Audible indicator beeps for a few seconds	Wait for unit to complete Power On Test, up to 30 seconds.
Start-up Period		The unit has recently been started and is in start-up mode, the output flow of the oxygen is temporarily < 82% Symbol will show yellow.	Yellow Startup Low O2 LED ON along with one O2% Status LED	No audible alarm during start-up period	Wait for unit to finish start-up period, up to 15 minutes
Normal Operation		Normal device operation (O2 >= 85.0%), no errors, startup period over Symbol will show green	Green Good O2 LED ON	OFF	Use device as desired
Check Cabinet Air Filter Notification		Notification – Check Cabinet Air Filter. LED turns ON every 7 days for 12 hours then restarts the timing cycle.	Yellow Check Cabinet Air Filter LED ON	OFF	Check Cabinet Filter and clean (wash/dry) or replace as necessary.

Check Compressor Intake Filter Notification		Notification – Check Compressor Intake Filter. LED turns ON every 133 days for 12 hours then restarts the timing cycle.	Yellow Check Compressor Intake Filter LED ON	OFF	Check Compressor Intake Filter and replace if necessary
Low Oxygen Output Concentration		Error – O2 below 85.0% ($\geq 60.0\% < 85.0\%$), startup period over. Symbol will show yellow.	Yellow Low O2 LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Continue to use device. Refer to Troubleshooting to try to resolve problem. If problem not resolved, contact your service technician for assistance.
Critical Low Oxygen Output Concentration		Error – O2 below 60.0% ($< 60.0\%$), startup period over. Symbol will show red.	Red Critical Low O2 LED ON. Red Service Required LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Refer to Troubleshooting to try to resolve problem. If problem not resolved, contact your service technician for assistance and to arrange for servicing of the unit.
Alarm Acknowledge Button Off		Alarm Acknowledge button not active (either no current alarm condition or audible alarm is actively sounding).	Alarm Acknowledge button LED OFF	If no current alarm, OFF. If alarm condition is active, audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Correct the Alarm Condition if audible alert is active (beeping).
Alarm Acknowledge Button On		Alarm Acknowledge button active (symbol is illuminated indicating active alarm silenced).	Alarm Acknowledge button LED ON	OFF	Correct the Alarm Condition to deactivate audible alert and return to normal operation.
Over Temperature		Warning – O2S gas over temperature or internal over temperature. Symbol will show red.	Over Temperature LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Refer to Troubleshooting to try to resolve problem. If problem not resolved, contact your Technician for assistance and to arrange for servicing of the unit.
Low Flow Rate		Warning – flow rate too low (no flow or obstruction). Flow rate below low flow threshold for 30 seconds. Symbol will show red.	Low Flow LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Refer to Troubleshooting to try to resolve problem. If problem not resolved, contact your Technician for assistance and to arrange for servicing of the unit.
High Flow Rate		Warning – flow rate too high. Error due to corrupted settings, OSD communication failure, fan not running, motor error, or rotary valve error. Symbol will show red.	Red Service Required LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Refer to Troubleshooting to try to resolve problem. If problem not resolved, contact your Technician for assistance and to arrange for servicing of the unit.
Device Malfunction		The device is experiencing a malfunction that requires servicing to correct. Symbol will show red.	Red Service Required LED light is illuminated	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Contact your service technician for assistance and to arrange for servicing of the unit
Power Failure		Power Failure Alarm – AC power not present, microcontroller running from a backup capacitor. Symbol will show red.	Red Service Required LED is flashing	The audible alarm is beeping intermittently. Alarm Acknowledge button is not available during a power failure alarm.	Turn unit off until electric power is restored.

4.1 INDICATORS AND ALARMS

Initial Inspection

1. Upon receiving, examine the unit for external damage. If the unit appears to have external damage, please contact Drive DeVilbiss for assistance.
2. Check to be sure the cabinet air filter and the compressor intake filter are in place.
3. Connect the power cord and plug the unit into an electrical outlet, turn the unit ON and check the audible/visual alarms. When the unit is turned ON, as part of the normal start-up process, all ten (10) LEDs on the front panel should illuminate and the audible alarm should sound when the unit is first turned ON. If ANY front panel LEDs DO NOT illuminate or the audible alarm DOES NOT sound, the alert system is not functioning properly. Refer to the Troubleshooting chart and contact your service technician if necessary.
4. Set the flow meter at the maximum recommended flow rate and allow the unit to run for 20 minutes. The internal oxygen sensor monitors the oxygen purity. If the oxygen is within specification, the Green Good Oxygen LED will be illuminated. If the Yellow Low Oxygen LED is illuminated, refer to Troubleshooting in the IFU, the Service Manual, or contact Drive DeVilbiss for assistance.

- With unit still running, unplug to test the power fail alarm. If the power fail alarm does not provide an audible alert, refer to the Service Manual or contact Drive DeVilbiss for assistance.

Preventive Maintenance Guidance

- Scheduled maintenance should be performed in accordance with the Preventive Maintenance/Service Guide table on page 16 or Between Patient Uses.
- Dispose of all accessories that are not suitable for reuse including oxygen tubing and connectors, cannula/mask, oxygen outlet connector, and humidifier bottle. Reusing accessories marked for single use may cause risk of cross contamination or loss of performance.
- Clean and inspect cabinet air filter and follow the Cleaning and Disinfection Instructions
- Clean the concentrator cabinet and inspect/replace filters in accordance with the table on page 16.
- Inspect all plugs, cords, and components. Replace any damaged or worn components.
- Check oxygen concentration with a calibrated oxygen analyzer and record the oxygen percentage. If the concentration is not within specification, refer to the troubleshooting section of the IFU or the Service Manual.
- Record the unit hours of use.
- Verify Audible Alert and Indicator LEDs at each service at startup and while operating.
- With unit still running, unplug to test the power fail alarm. If the power fail alarm does not provide an audible alert, refer to the Service Manual.**

NOTE

Scheduled maintenance should be performed in accordance with the Preventive Maintenance/Service Guide table above or Between Patient Uses.

4.2 TROUBLESHOOTING

TROUBLESHOOTING GUIDE

SYMPTOM	POSSIBLE CAUSE	REMEDY
A. Unit does not operate. All lights are off when the power switch is On. Audible alert is pulsing	1. Power cord not properly inserted into wall outlet.	1. Check power cord connection at the wall outlet. Also check the mains connection on the back of the unit.
	2. No power at wall outlet.	2. Check your circuit breaker and reset if necessary. Use a different wall outlet if the situation occurs again. If the above remedies do not work, contact your service technician.
B. Yellow Start Up Low O2 LED is on.	1. Unit in "start up" mode.	1. Allow unit up to 15 minutes to complete start up period.
C. Yellow Low Oxygen LED is on and intermittent audible signal is sounding.	1. Flow meter is not properly set.	1. Ensure the flow meter is properly set to the prescribed number.
	2. Filter door vents are blocked.	2. Check filter door vents and ensure that the openings are not blocked.
	3. Exhaust is blocked.	3. Check the exhaust area and make sure there is nothing restricting the unit exhaust. If the above remedies do not work, contact your service technician.
D. Red Critical Low O2 LED and Red Service Required LEDs are illuminated and an intermittent audible signal is sounding.	1. Flow meter is not properly set.	1. Ensure the flow meter is properly set to the prescribed number.
	2. Filter door vents are blocked.	2. Check filter door vents and ensure that the openings are not blocked.
	3. Exhaust is blocked.	3. Check the exhaust area and make sure there is nothing restricting the unit exhaust. If the above remedies do not work, contact your service technician.
E. Unit operates. Red Over Temperature LED is illuminated. Audible alert may be sounding.	1. Filter door vents are blocked.	1. Check filter door vents and ensure that the openings are not blocked.
	2. Exhaust is blocked.	2. Check the exhaust area and make sure there is nothing restricting the unit exhaust. If the above remedies do not work, contact your service technician.

SYMPTOM	POSSIBLE CAUSE	REMEDY
F. Unit operates. Red Low Flow LED is illuminated. Audible alert may be sounding.	1. Blocked or defective cannula, face mask, or oxygen tubing.	1. Detach cannula or face mask. If proper flow is restored, clean or replace if necessary. Disconnect the oxygen tubing at the oxygen outlet. If proper flow is restored, check oxygen tubing for obstructions or kinks. Replace if necessary.
	2. Blocked or defective humidifier bottle.	2. Detach the humidifier from the oxygen outlet. If proper flow is obtained, clean or replace humidifier.
	3. Flow meter set too low.	3. Set flow meter to prescribed flow rate. If the above remedies do not work, contact your service technician.
G. Red Service Required LED is illuminated and an intermittent audible signal is sounding.	1. Electronic assembly or other internal component malfunction.	1. Turn your unit OFF. Call your service technician and switch to alternate oxygen source
H. Unit operates. Audible low-frequency vibration sound is detected.	1. Electronic assembly malfunction.	1. Turn your unit OFF. Call your service technician and switch to alternate oxygen source.
I. Unit operates. Any of the visual and audible alerts do not function when the device is turned ON	1. Electronic assembly malfunction.	1. Turn your unit OFF. Call your service technician and switch to alternate oxygen source.
J. Red Service Required LED flashes, no AC power present, with audible alarm	1. Power Failure Alarm	1. Turn unit off until electricity or power is restored.
K. If any other problems occur with your oxygen concentrator.		1. Turn your unit OFF. Call your service technician and switch to alternate oxygen source.

5.0 Specifications

PulmO2 10 LITER OXYGEN CONCENTRATOR	
Catalog Number	1060AW
Delivery Rate	1 to 10 LPM
Maximum Recommended Flow (@ nominal outlet pressures of zero and 7 kPa)**	10 LPM
Outlet Pressure	58.6 ± 3.5 kPa (8.5 ± 1.0 PSIG)
Oxygen Percentage	1-10 LPM = 93%-6%, +3%
Minimum volume fraction of delivered oxygen at the maximum rated flowrate	85% O ₂ @ 10 LPM
Flow Meter Accuracy	+/- 0.5 SLPM, at 10 SLPM; +0, -0.5 SLPM
5.1 ENVIROMENTAL	
Operating Atmospheric Pressure	1013 hPa to 795 hPa 0-2000 meters (0-6562 feet) Tested at nominal voltage.
Operating Environment Range	5°C (41°F) – 40°C (104°F)
Operating Relative Humidity Range	15% to 95%, non-condensing
Power Consumption	230 VAC, 50Hz: 520 Watts Maximum
Sound Pressure Level at 3 and 10 LPM	<60 dBA
Sound Power Level at 3 and 10 LPM	<70 dBA
5.2 PHYSICAL	
Weight	25.2 kg (55.5 lbs)
Dimensions	70.5 cm H x 40.6 cm W x 35.6 cm D (27.75" H x 16" W x 14" D)
Transport and Storage Conditions	-25°C - (-13°F) to 70°C (158°F), humidity range of 15% to 95% non-condensing

PulmO2 10 LITER OXYGEN CONCENTRATOR

Catalog Number	1060AW
5.3 STANDARDS COMPLIANCE	
Approval Body and Safety Standard	TUV IEC 60601-1:2020 IEC 60601-1-6:2020 IEC 60601-1-11:2020 EN ISO 80601-2-69:2020
EMC Compliance To	EN60601-1-2:2020
5.4 ELECTRICAL	
Electrical Rating	100-240 VAC, 50 / 60Hz, 3.2 Amp
5.5 PRESSURE	
Maximum Limited Pressure under Normal condition under Single Fault condition	186 kPa (27 PSIG) 276 kPa (40 PSIG)
Operating System	Time Cycle / Pressure Swing
5.6 CLASSIFICATION	
Equipment Class and Type	Class II Equipment; Type BF Applied Part IP21 Ingress Protection

Specifications subject to change without notice.

Units of measure for volumes, flows and leakages are expressed in STPD: 101,3 kPa at an operating temperature of 20 °C, dry

Oxygen Concentration by Flow Rate

Flow L/m	%O ₂
10	93% -6%, +3%
9	93% -6%, +3%
8	93% -6%, +3%
7	93% -6%, +3%
6	93% -6%, +3%
5	93% -6%, +3%
4	93% -6%, +3%
3	93% -6%, +3%
2	93% -6%, +3%
1	93% -6%, +3%

5.7 Service Life

The expected service life of the PulmO2 10 L oxygen concentrator, which includes the performance of any required service is 5 years. Drive DeVilbiss recommends for optimal service life that the PulmO2 10L Oxygen Concentrator be operated for at least 30 minutes after it is powered on. Shorter periods of operation, operating in extreme temperature/humidity conditions or in the presence of contaminants, and/or handling and storage conditions outside those specified, may affect the long term reliable operation of the product.

5.8 Return and Disposal

This device may not be disposed of with household waste. After use of the device, please return the device to the provider for disposal. This device contains electrical and/or electronic components that must be recycled per EU Directive 2012/19/EU-Waste Electrical and Electronic Equipment (WEEE).

Store used sieve beds in sealed airtight container to minimize release of absorbed substances. Waste must be handled in accordance with all federal, state, provincial, and local regulations. There are no known hazardous contaminants absorbed into the molecular sieve therefore used sieve beds may be disposed of as residential waste. Never dispose by means of public sewers or drainage.

The disposal of infectious accessories (e.g. nasal cannula from an infected user) must be made via an approved waste disposal company. Names and addresses can be obtained from the local municipality.

IMPORTANT

This device contains electrical and/or electronic equipment. Follow local governing ordinances and recycling plans regarding disposal of device components.

APPENDIX A

Electromagnetic Compatibility Information

IMMUNITY TEST LEVELS			
TEST DESCRIPTION	SPECIFICATION	NOTES	RESULTS
CISPR11 & IEC 60601-1-2:2014 + AMD1:2020 (EDITION 4.1)			
RF Emissions	CISPR 11	Group 1	Complies
RF Emissions	CISPR 11	Class B	Complies
Harmonic Current Emissions	EN 61000-3-2:2014	AC Input Class B	Complies
Voltage Fluctuations & Flicker Emissions	EN 61000-3-3:2013	AC Input	Complies
Electrostatic Discharge Immunity (ESD)	EN 61000-4-2:2008	±15kV (Air) ±8kV (Contact)	Complies
Radiated Electromagnetic Field Immunity	EN 61000-4-3:2006	80MHz to 2.7GHz @ 10V/m 2.7-6.0GHz @ 3V/m 80% AM @ 1kHz	Complies
Fast Transient/Burst Immunity	EN 61000-4-4:2004	±2kV for Power Lines ±1kV for I/O Lines	Complies
Surge Immunity	EN 61000-4-5:2006	±1kV differential ±2kV common	Complies
Conducted RF Immunity	EN 61000-4-6:2009	3Vrms 150kHz to 80MHz 6Vrms ISM Band frequencies between 150kHz-80MHz	Complies

Maximum length of AC line cord = 3m (9.8 ft)

WARNINGS

- ▶ Do not bring the device or accessories into a Magnetic Resonance (MR) environment as it may cause unacceptable risk to the patient or damage to the oxygen concentrator or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.
- ▶ Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the oxygen concentrator. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, disconnect the power cord and discontinue use. Contact your healthcare provider or clinical staff.
- ▶ This device is suitable for use in home and healthcare environments except for near active HF surgical equipment and the RF shielded room of a magnetic resonance imaging system, where the intensity of electromagnetic disturbances is high.
- ▶ Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- ▶ Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the oxygen concentrator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- ▶ Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

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DEVILBISS® AND OSD® ARE REGISTERED TRADEMARKS OF DEVILBISS HEALTHCARE.
© 2024 DeVilbiss Healthcare LLC. 10.24 All Rights Reserved.
Drive DeVilbiss Healthcare Ltd.
Drive DeVilbiss Australia Pty Limited
DeVilbiss Healthcare GmbH
DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501 • USA