DRIVE DEVILBISS® PULMO2 10-LITER OXYGEN CONCENTRATOR SERVICE MANUAL

DANGER – NO SMOKING

MODEL 1060AW

R_x ONLY MD

Device is not manufactured with natural rubber latex.

Assembled in USA

TABLE OF CONTENTS

GENERAL INFORMATION

Introduction.....	3
Symbol Definitions.....	3
Important Safeguards.....	4

UNPACKING AND SETUP

Initial Inspection.....	8
Patient Setup.....	8
Operating Instructions.....	8

MAINTENANCE

Patient Alert System.....	9
Alarm Function Testing.....	10
Service Life.....	11
Routine Patient Maintenance.....	11
Periodic Homecare Provider Preventative Maintenance.....	11
Preventative Maintenance Summary.....	12
Provider's Notes.....	13
Service and Maintenance Guidance.....	13

TROUBLESHOOTING

System Operation.....	14
Normal Operating Sequence.....	14
Control Panel Alarms Overview.....	15
Overview Of Alarms And Service Indicators.....	15
Simplified Troubleshooting.....	16
Troubleshooting Chart.....	17

COMPONENT TESTING, REPAIR, AND REPLACEMENT

Proper Repair Procedures.....	18
Cabinet Removal.....	18
Accumulator Tank.....	19
Accumulator Pressure Test.....	19
Accumulator Tank Removal.....	19
Check Valves.....	20
Compressor.....	20
Cooling Fan.....	23
Flow Meter.....	24
Fused IEC Connector.....	24
Hour Meter.....	24
Leak Testing.....	24
Manifold.....	24
Molecular Sieve Beds.....	24
Power Cord.....	25
Power Switch.....	25
Pressure Regulator.....	25
Printed Circuit Board (PC Board).....	26
Rotary Valve.....	26

FIGURES, DIAGRAMS AND PARTS LIST

1060AW Unit.....	28
Accessories and Tools.....	29
Pneumatic and Wiring Diagrams.....	31

ORDERING INFORMATION AND PARTS RETURN

Ordering Information.....	33
Returns.....	33

WARRANTY

Return and Disposal.....	34
--------------------------	----

SPECIFICATIONS

ELECTROMAGNETIC COMPATIBILITY INFORMATION.....	36
--	----

GENERAL INFORMATION

INTRODUCTION

This service manual was designed to provide Biomedical technicians and clinicians with the proper maintenance, service, safety, and repair procedures for the Drive DeVilbiss Oxygen Concentrator.

Read and understand all the information contained in this service manual before attempting to operate or perform any maintenance on the concentrator.



An oxygen concentrator is a device that delivers highly concentrated oxygen for therapeutic applications.

Room air is a mixture of 78% nitrogen, 21% oxygen, 1% argon and other gases. The concentrator draws in room air, separates the nitrogen from the oxygen, and delivers concentrated oxygen to the patient through an oxygen port.

NOTE– DeVilbiss reserves the right to alter or change the design of the Drive DeVilbiss Oxygen Concentrator series. Hence, slight differences in construction or components may exist between the unit in hand and what is described in this manual.

SYMBOL DEFINITIONS

	Manufacturer		Medical Device		Serial Number		LOT Number		Clinical Staff
	Catalog Number		Class II Equipment		Attention Required		Recurring Reminders		Service Technician
	Humidifier Bottle		Place on a hard surface		Well ventilated area		30 cm away from wall		On Off
	Operating Temperature Range +5 to +40°C (+41 to +104°F)		Operating Atmospheric Pressure Range 795 to 1013 hPa (Sea level to 2000m (6562 ft))		TUV Rheinland C-US approval mark		Type BF applied part		User
	Storage Humidity Limit		Keep Dry		Storage Temperature Limits		Hour Meter		Alternating Current
	MRI Unsafe symbol – Unsafe for Magnetic Resonance Environment <i>This symbol has a red circle and diagonal bar on the product label.</i>		Do not use wet cabinet air filter <i>This symbol has a red circle and diagonal bar on the product label.</i>		Do not cover <i>This symbol has a red circle and diagonal bar on the product label.</i>		Do not Disassemble <i>This symbol has a red circle and diagonal bar on the product label.</i>		HF RF
	CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.		Over Temperature LED <i>This symbol illuminates Red when active</i>		Startup Low O2 Indicator, <i>Oxygen is low during initial start up.</i> <i>This symbol illuminates Yellow when active.</i>		Alarm Acknowledge Button and Active LED <i>This symbol illuminates White when active.</i>	DO NOT USE NEAR SURGICAL EQUIPMENT with High Frequency Surgical Equipment / Radio Frequency Equipment <i>This symbol has a red circle and diagonal bar on the product label.</i>	
	Check Cabinet Air Filter Notification LED <i>This symbol illuminates Yellow when active</i>		Check Compressor Intake Filter Notification LED <i>This symbol illuminates Yellow when active</i>		Low Flow LED <i>This symbol illuminates Red when active</i>		Service Required LED <i>This symbol illuminates Red when active.</i>		Oxygen Outlet marked with rated range of gas pressure and rated range of gas flowrate

	It is mandatory to read and understand the operating instructions prior to use. <i>This symbol has a blue background on the product label.</i>		Electric Shock Hazard. Cabinet to be removed by authorized personnel only. <i>This symbol has a yellow background on the product label.</i>		No smoking near patient or device. <i>This symbol has a red circle and diagonal bar on the product label.</i>		Use no Oil, Grease or Lubricants. <i>This symbol has a red circle and diagonal bar on the product label.</i>		Do not use near heat or open flames <i>This symbol has a red circle and diagonal bar on the product label.</i>
	This device contain electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU- Waste Electrical and Electronic Equipment (WEEE)		General Warning <i>This symbol is used throughout this manual to indicate hazardous situations to avoid.</i>		Important Information <i>This symbol is used throughout this manual to indicate important information you should know.</i>		Note and Information Symbol <i>This symbol is used throughout this manual to indicate notes, useful tips, recommendations and information.</i>		Ingress Protection - Protected against finger access to hazardous parts; protected against vertically falling water drops.
 O2% (Oxygen Percentage) Status LEDs	 Maximum recommended flow rate: 10 LPM		QR code link to user and service information		Caution is necessary when using the device.		Consult instructions for use		
	Unique Device Identifier								

IMPORTANT SAFEGUARDS

Read this entire guide before using your Drive DeVilbiss concentrator. Important safeguards are indicated throughout this guide. Pay special attention to all safety information. Imminently and potentially hazardous information is highlighted by these terms:

DANGER

Indicates an imminently hazardous situation which could result in death or serious injury to the user or operator if not avoided.

WARNING

Indicates a potentially hazardous situation which could result in death or serious injury to the user or operator if not avoided.

CAUTION

Indicates a potentially hazardous situation which could result in property damage, injury, or device damage if not avoided.

IMPORTANT

Indicates important information you should know.

NOTE

Indicates notes, useful tips, recommendations, and information.

READ ALL INSTRUCTIONS BEFORE USING.

GENERAL INFORMATION

DANGER

1. NO SMOKING signs should be prominently displayed.
2. Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy.
3. Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located.
 - If you intend to smoke, you must always turn the oxygen concentrator off, remove the cannula and leave the room where either the cannula or mask or the oxygen concentrator is located. If unable to leave the room, you must wait 10 minutes after you have turned off the oxygen concentrator before smoking.
4. Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions if the oxygen concentrator is turned on but not in use. The oxygen will make the materials flammable. Turn the oxygen concentrator off when not in use to prevent oxygen enrichment.
5. Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.
6. Open flames during oxygen therapy are dangerous and are likely to result in fire or death. Do not allow open flames within 2 m (6.5 feet) of the oxygen concentrator or any oxygen carrying accessories.
7. Drive DeVilbiss oxygen concentrators are equipped with a fire mitigating outlet fitting that prevents propagation of fire into the unit.

WARNING

1. To avoid electric shock, do not plug the concentrator into an AC outlet if the concentrator cabinet is broken. Do not remove the concentrator cabinet. The cabinet should only be removed by a qualified Drive DeVilbiss technician. Do not apply liquid directly to the cabinet or utilize any petroleum-based solvents or cleaning agents.
2. Improper use of the power cord and plugs can cause a burn, fire or other electric shock hazards. Do not use the unit if the power cord is damaged.
3. Ensure the mains power cord is fully inserted into the concentrator connector and the power cord plug is completely inserted into a fully functioning AC wall outlet. Failure to do so may cause an electrical safety hazard.
4. The accessories (nasal cannula, masks, oxygen tubing, humidifiers, etc.) that supply oxygen to the patient must be equipped with a means that, in case of fire, stops the propagation of fire through the accessory for the safety of the patient and others. A fire activated flow-stop or thermal fuse device, if available, should be used with the oxygen supply accessories. These types of flow-stop devices stop the flow of oxygen to the patient in the event of fire. This means of fire protection should be located as close to the patient as practicable.
5. Locate oxygen tubing and power supply cords to prevent tripping hazards and reduce the possibility of entanglement or strangulation.
6. Do not lubricate fittings, connections, tubing or other accessories of the oxygen concentrator to avoid the risk of fire and burns.
7. Do NOT use lubricants, oils or grease.
8. Before attempting any cleaning procedures, turn the unit off.
9. Use only water-based lotions or salves that are oxygen-compatible before and during oxygen therapy. Never use petroleum or oil-based lotions or salves to avoid the risk of fire and burns.
10. Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.
11. If you feel discomfort or are experiencing a medical emergency while undergoing oxygen therapy, seek medical assistance immediately to avoid harm.
12. Geriatric, pediatric or any other patient unable to communicate discomfort can require additional monitoring and/or a distributed alarm system to convey the information about the discomfort and/or the medical urgency to the responsible caregiver to avoid harm.
13. Use of this device at an altitude above 10,000 feet (3050 meters) or above a temperature of 104°F (40°C) or greater than 95% relative humidity may affect the flow rate and the percentage of oxygen and consequently the quality of the therapy. Refer to specifications for details regarding parameters tested.
14. To ensure you receive the therapeutic amount of oxygen delivery according to your medical condition, the Oxygen Concentrator must:
 - be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels.
 - be used with the specific combination of parts and accessories that are in line with the specification of the concentrator manufacturer and that were used while your settings were determined.
15. Your delivery settings of the oxygen concentrator should be periodically reassessed for the effectiveness of therapy.
16. For your safety, the oxygen concentrator must be used according to the prescription determined by your physician.
17. Under certain circumstances, oxygen therapy can be hazardous. Seek medical advice before using an oxygen concentrator.
18. The proper placement and positioning of the prongs of the nasal cannula in the nose is critical to the amount of oxygen delivered to the respiratory system of the patient.
19. Do not modify the oxygen concentrator or accessories.
20. Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
21. This device may be affected by children unintentionally changing the device setting or inserting foreign objects into device openings. Supervise children when the device is in operation.
22. The device can be affected by pet hair or insects getting into the device or obstructing cooling vents on the device.
23. Concentrated oxygen flow rates above 1 LPM for Neonatal and 2 LPM for infant patients can create a risk of high FiO₂ (>50%) which may cause issues like retinopathy of prematurity and bronchopulmonary dysplasia. In infants, very high FiO₂ may lead to oxygen toxicity - damage to lung tissue.

GENERAL INFORMATION

WARNING

MR Unsafe

1. Do not bring the device or accessories into a Magnetic Resonance (MR) environment as it may cause unacceptable risk to the patient or damage to the oxygen concentrator or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.
2. Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the oxygen concentrator. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, disconnect the power cord and discontinue use. Contact your home care provider.
3. This device is suitable for use in home and healthcare environments except for near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of Electromagnetic DISTURBANCES is high.
4. Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
5. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the oxygen concentrator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
6. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
7. Do not use ACCESSORIES, detachable parts and materials not described in the instructions for use.
8. Do not interconnect this oxygen concentrator with other equipment not described in the instructions for use.

CAUTION

Risk of injury or damage

- Operate the unit in a cool, dry area with good ventilation, located on a hard surface, avoid thick rugs or carpeting. NEVER block the air intake or exhaust vents.
- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Use only recommended cleaning and disinfection agents to avoid damage to the device.
- Keep exposed body parts, such as hands and feet, a minimum of 46 inches (1.2 meters) away from the exhaust vents to avoid the risk of burns. Single fault conditions may result in visual and audible alerts and alarms. The oxygen concentrator should be located in a well-ventilated area. DO NOT operate the unit in a closed or confined space, such as a closet, bathroom, etc. Avoid operating the device near smoke pollutants and fumes.
- To prevent product damage, do not attempt to operate the unit without the air filter or while the filter is still damp.
- Incompatible parts or accessories can result in degraded performance. Use only recommended parts and filters to ensure reliable operation of the product.
- To prevent product damage, do not attempt to operate the unit without the compressor intake filter.
- Use of harsh chemicals (including alcohol) is not recommended. If bactericidal cleaning is required, a non-alcohol based product should be used to avoid inadvertent damage
- Install, use and maintain the PulmO2 Series 10 Liter Oxygen Concentrator following this instruction guide to minimize the environmental impact of the oxygen concentrator through its expected life. The oxygen concentrator consumes maximum 520 watts during normal use. Water can be used as needed for humidification of the delivered oxygen. The Compressor Intake Filter and Final Hepa Filter are disposable and could need to be replaced twice during the expected service life. Other byproducts created during the normal use of the oxygen concentrator are heat, noise and nitrogen gas. Follow instructions to minimize the effects of heat and noise. Nitrogen gas discharged by the oxygen concentrator is quickly dissipated back into room air. The oxygen concentrator contains no hazardous substances.

Table 1 - Model 1060AW

Description	Maximum Temperature		Max safe contact time
	1060AW Series		
	°F	°C	
Outlet Vents	145.8	63.2	Less than 1 minute
Oxygen outlet fitting	107.6	42.0	Less than 10 minutes
Oxygen Tubing Coupler	106.0	41.1	Less than 10 minutes
Cannula tubing at the outlet fitting	106.2	41.2	Less than 10 minutes

CAUTION

To avoid possible injury or loss of therapy, turn the device off and disconnect from patient before moving the device.

GENERAL INFORMATION

IMPORTANT

- It is recommended that the clinical staff or provider lock the flow control knob to prevent inadvertent adjustment. A flow setting other than prescribed may affect the patient therapy.
- Do not service or clean this device while in use with a User.
- The 1060AW Series oxygen concentrator, its parts and accessories are specified for use at flow rates between 1 LPM and 10 LPM.
- When used in a homecare environment, the patient is the intended operator.

SAVE THESE INSTRUCTIONS.

UNPACKING AND SETUP

INITIAL INSPECTION

1. Upon receiving, examine the unit for external damage. If the unit appears to have external damage, please contact DeVilbiss for assistance.
2. Check to be sure the cabinet air filter and the compressor intake filter are in place.
3. Plug the unit into an electrical outlet, turn the unit ON and check the audible/visual alarms. When the unit is turned ON, as part of the normal start-up process, all ten (10) LEDs on the front panel should illuminate and the audible alarm should sound when the unit is first turned ON. If ANY front panel LEDs DO NOT illuminate or the audible alarm DOES NOT sound, the alert system is not functioning properly. Refer to the Troubleshooting chart and contact your biomedical technician if necessary.
4. Set the flow meter at the maximum recommended flow rate and allow the unit to run for 20 minutes. The internal oxygen sensor monitors the oxygen purity. If the oxygen is within specification, the Green Good Oxygen LED will be illuminated. If the Yellow Low Oxygen LED is illuminated, refer to Troubleshooting in the IFU, the Service Manual, or contact DeVilbiss for assistance.
5. With unit still running, unplug to test the power fail alarm. If the power fail alarm does not provide an audible alert, refer to the Service Manual or contact DeVilbiss for assistance.

When ready for operation

WARNING

The accessories (nasal cannula, masks, oxygen tubing, humidifiers, etc.) that supply oxygen to the patient must be equipped with a means that, in case of fire, stops the propagation of fire through the accessory for the safety of the patient and others. A fire activated flow-stop or thermal fuse device, if available, should be used with the oxygen supply accessories. These types of flow-stop devices stop the flow of oxygen to the patient in the event of fire. This means of fire protection should be located as close to the patient as practicable.

1. Attach the nasal cannula to the oxygen tubing per the manufacturer's directions
2. Follow the Operating Instructions.

OPERATING INSTRUCTIONS

1. Remove the power cord completely from the strap. Make sure the power switch is in the "Off" position.
2. Ensure cord is connected to the unit before inserting plug into an appropriate electrical outlet.

WARNING

Improper use of the power cord and plugs can cause a burn, fire, or other electric shock hazards. Do not use the unit if the power cord is damaged.

Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy. Keep the oxygen concentrator and cannula at least 6.5 feet (2 m) from hot, sparking objects or naked sources of flame.

WARNING

The Drive DeVilbiss 1060AW oxygen concentrator is equipped with a fire mitigating outlet fitting that prevents propagation of fire into the unit.

The accessories (nasal cannula, masks, oxygen tubing, humidifiers, etc.) that supply oxygen to the patient must be equipped with a means that, in case of fire, stops the propagation of fire through the accessory for the safety of the patient and others. A fire activated flow-stop or thermal fuse device, if available, should be used with the oxygen supply accessories. These types of flow-stop devices stop the flow of oxygen to the patient in the event of fire. This means of fire protection should be located as close to the patient as practicable.

WARNING

As part of the normal start-up process, all ten (10) LEDs on the front panel should illuminate and the audible alarm should sound when the unit is first turned ON. If ANY front panel LEDs DO NOT illuminate or the audible alarm DOES NOT sound, the alert system is not functioning properly. Refer to the Troubleshooting chart and contact your biomedical technician if necessary.

3. Press the power switch to the "ON" position. When the unit is turned ON, all ten (10) LEDs on the front panel will illuminate briefly and an audible signal will briefly alarm confirming that the LEDs and audible signal are functioning properly. The unit will then operate in "start-up" mode with the Startup Low O₂ LED illuminated until a normal oxygen level is achieved, at which time the Green Good O₂ LED illuminates and will remain lit. The start-up mode may take up to 15 minutes..

NOTE— Drive DeVilbiss recommends for optimal service life that the Drive DeVilbiss Oxygen Concentrator to be operated for at least 30 minutes after it is powered on. Shorter periods of operation, operating in extreme temperature/humidity conditions or in the presence of contaminants, and/or handling and storage conditions outside those specified, may affect the long term reliable operation of the product.

4. Slowly turn the flow meter knob until the flow meter ball is centered on the line next to the appropriate flow rate.

NOTE— When the flow meter knob is turned clockwise, the flow decreases (and eventually will shut off the oxygen flow). When the knob is turned counter-clockwise, the flow increases.

NOTE— For prescriptions of 10 LPM, be sure the ball is centered on the 10 liter line; the ball should not touch the red line. Setting the flow higher than 10 may cause the oxygen purity level to drop.

NOTE— The low-flow alarm may activate if the flow meter ball is set below 1 lpm. The unit will continue to run; however, the Service Required light will come on accompanied by an audible alarm. Adjust the flow meter to your prescribed flow.

NOTE— The unit may require up to 20 minutes for the oxygen concentration and flow rate to stabilize. The flow rate should be monitored and readjusted if necessary.

5. The flow meter has a locking device. If it is necessary to preset and lock in the prescribed flow rate, tighten the set screw located on the hex nut just below the control knob using a 1/16" Allen bit. No adjustment can be made without loosening the set screw.
6. The Drive DeVilbiss oxygen concentrator is now ready for use.

PATIENT ALARM SYSTEM

The Drive DeVilbiss Oxygen Concentrator patient alarm system will detect unit component failure. This system is comprised of both visible and audible alarms which signal the patient if a malfunction should occur.

DeVilbiss OSD® Operation

The OSD (Oxygen Sensing Device) is a device within Drive DeVilbiss concentrators that monitors the oxygen produced by the unit. It also monitors the flow.

The OSD lights on the top panel are defined as follows:

Yellow Start-Up Low Oxygen LED - the unit has just been started and hasn't yet reached the acceptable oxygen level.

Green Good Oxygen LED – acceptable oxygen level.

Yellow Low Oxygen LED – below the good oxygen level.

Red Critical Low Oxygen LED - below the critical oxygen level.

Low Flow Rate LED - flow rate is too low.

If the oxygen purity falls below the acceptable level, the green Good Oxygen LED will shut off, the yellow Low Oxygen LED will illuminate, and an intermittent audible signal will sound.

If the oxygen purity continues to fall into the critical low level, the Red Critical Low Oxygen LED will illuminate, along with an intermittent audible alarm.

If the flow rate is too low or there is no flow, the Low Flow Rate LED will illuminate, and an intermittent audible signal will sound.

Refer to the Troubleshooting section in this guide. Do not attempt any other maintenance. Contact your biomedical technician immediately if the issue is not resolved.

NOTE– Refer to the Alerts section below for specific alarm settings.

NOTE– After power on, the electronics continuously monitor the oxygen sensor. If a fault is detected, the appropriate LED(s) and/or audible alarms will activate.

The first 15 minutes, the unit will be in Start Up mode. The oxygen purity is continuously monitored and the yellow Start Up Low Oxygen LED will be On. The green Good Oxygen LED will turn on as soon as the therapeutic oxygen levels are obtained. After 15 minutes stabilization time, if the O₂ is less than 82% the yellow Low O₂ LED will illuminate and a beeping audible alarm will occur. If the oxygen level falls below 60% (after startup) then the red Critical Low Oxygen LED will illuminate along with a beeping audible alarm.

Alarms:

There are several visible service alarms located on the control panel. The audible alarm system is internally powered; no batteries are required. If the indicator lights illuminate or the audible alarm sounds other than during start-up, a problem has occurred. Detailed Alarm Information is found in Troubleshooting section of this manual.

- Power Failure (Red Service Required LED and intermittent audible alarm)
- Low Flow Rate (Low Flow LED and intermittent audible alert)
- Over Temperature (Over Temperature LED and intermittent audible alarm)

Below normal O₂:

- The yellow Low Oxygen LED will illuminate with an intermittent audible alarm at approximately <82%.
- The red Critical Low Oxygen LED will illuminate with an intermittent audible alarm at <60%.

The visible and audible alarms will activate for a minimum of two minutes in a no power situation. If the unit is turned On without power or power is removed later, the alarm will sound within 10 seconds. After that time, the alarm will produce an audible pulse every few seconds while the LED flashes. Power for this alert is provided by a capacitor on the PC board.

NOTE– If the concentrator has been unused for an extended period, the unit must run for a minimum of 2 minutes before the power fail alarm will activate. 10 minutes of run time is needed to fully charge the capacitors to run the power fail alarm for the entire 120 second specification.

The PC (printed circuit) board is responsible for controlling the system and alarms.

NOTE– A high pressure condition is indicated by the audible (a “popping” sound) release of pressure from a pressure relief valve located on the compressor head.

NOTE– Settings below 1 LPM may activate the low flow alarm.

ALARM FUNCTION TESTING

The 1060AW series is designed to activate alarms when certain conditions or failures occur. The alarm functions may be tested following the procedures below:

1. **Overheating:**
 - a. Remove the front and rear covers from the concentrator; then disconnect the cooling fan from the printed circuit board.
 - b. Replace the front and rear covers.
 - c. Place the concentrator in a location that has an ambient temperature of approximately 70° F. (21°C) Then plug the unit into the appropriate mains voltage and turn it on.
 - d. Allow the unit to operate until the Over Temperature Alarm activates, which should be within approximately two hours.
2. **Compressor Failure:**
 - a. Remove the rear cover from the concentrator; then disconnect the compressor electrical connector from the main wire harness.
 - b. Plug the unit into the appropriate Mains voltage and turn it on.
 - c. Allow the unit to operate until the Motor Error activates, which should be within approximately two minutes.
3. **Low Flow / Obstruction of Gas Pathway:**
 - a. Plug the concentrator into the appropriate mains voltage and turn it on.
 - b. Allow the device to run for several minutes.
 - c. Turn the flow meter off so that there is no oxygen flowing out of the unit.
 - d. Allow the unit to operate until the Low Flow Rate Alarm activates.
 - e. Increase the flow to 2 LPM and confirm that the alarm condition ends.
4. **Oxygen Generation Mains Failure:**
 - a. Remove the front and rear covers from the concentrator.
 - b. Disconnect the tubing from the bottom of the flow meter.
 - c. Plug the unit into the appropriate mains voltage and turn it on.
 - d. Turn the flow meter to 10 LPM.
 - e. Allow the unit to operate until the Service Required Alarm activates.
5. **Power Supply Failure**
 - a. Connect the oxygen concentrator to AC power and turn the power switch on. Allow the device to run for a minimum of 10 minutes.
 - b. With the power switch in the On position, unplug the AC power cord from the outlet.
 - c. The alarm condition (blinking service required LED and intermittent audible alert) should occur and continue for a minimum of 120 seconds.
 - d. Reconnect the AC power cord and confirm that the alarm condition ends.
6. **Malfunction – O₂S Gas Temperature High**
 - a. This alarm condition is tested automatically during start-up.

7. **Malfunction – Corrupted Settings**
 - a. This alarm condition is tested automatically during start-up.
10. **Malfunction – Non-Recoverable Valve Error**
 - a. This alarm condition is tested automatically during start-up.
11. **Malfunction – O2S Oxygen Sensor Communication Failure**
 - a. This alarm condition is tested automatically during start-up.
12. **Low Oxygen Concentration – Startup Period Over**
 - a. This alarm condition is tested automatically during start-up. Once the oxygen concentration reaches 85% the green light comes on.

SERVICE LIFE

The expected service life of the 1060AW series oxygen concentrator, which includes the performance of any required service or maintenance, is 5 years. The expected service life is based on the operation of the device in accordance with all manufacturer guidance for safe use, maintenance, servicing, storage, shipping, handling, and general operation.

The actual service life of the unit, and in particular the service life of certain subcomponents, including the Filters, Sieve Beds and Compressor Cup Seals, will vary based on a number of variables, including the operating environment, storage environment, shipping, handling, performance of preventive maintenance, and both the frequency and intensity of use.

The 1060AW series oxygen concentrators have internal sensors and diagnostic systems designed to monitor the system performance, including the oxygen concentration (purity), flow and temperature. The 1060AW concentrators will alert the user when the device requires maintenance or service. Please see the Troubleshooting and Maintenance Sections for more detailed information.

ROUTINE PATIENT MAINTENANCE

CAUTION

Incompatible parts or accessories can result in degraded performance. Use only recommended parts and filters to ensure reliable operation of the product.

The oxygen patient should perform the following maintenance:

Cannula/Mask, Tubing, & Humidifier Bottle

The patient should clean and replace the cannula, tubing, and humidifier bottle according to the manufacturer's instructions.

Oxygen Outlet Connector

Remove the oxygen outlet connector and clean as needed or replace between patients.

Cabinet Air Filter

The gross particle air filter should be inspected when the Check Cabinet Filter LED is illuminated and cleaned as needed by the user or caregiver. Replace if torn or damaged. To clean, follow these steps:

NOTE– Frequency of inspection and cleaning of filter may be dependent upon environmental conditions like dust and lint.

1. Remove the air filter located in the door on the back of the unit.
2. Wash in a solution of warm water and dishwashing detergent.
3. Rinse thoroughly with warm tap water and towel dry. The filter should be completely dry before reinstalling.

CAUTION

To prevent product damage, do not attempt to operate the unit without the air filter or while the filter is still damp

NOTE– The air filter should be monitored more closely in environments with abnormal amounts of dust and lint.

CAUTION

Operation of the Drive DeVilbiss Oxygen Concentrator in extreme environments or without the cabinet air filter will prematurely occlude the compressor intake filter and cause a decrease in the unit performance.

Compressor Intake Filter

The compressor intake filter should be inspected when the Check Compressor Intake Filter LED illuminates and replaced as necessary. To check or replace, follow these steps:

1. Open the filter door on the rear of the concentrator.
2. Remove the filter from behind the door and inspect the color and debris. If filter is discolored to a dark gray, it should be replaced.

CAUTION

To prevent product damage, do not attempt to operate the unit without the compressor intake filter.

Exterior Cabinet

The patient should clean the concentrator exterior cabinet weekly by using a damp cloth or sponge with a mild household cleaner and wiping it dry. Inspect the filter door vents periodically, and wipe with a dry cloth as needed to remove dust.

WARNING

Do not apply liquids directly to the cabinet or utilize any petroleum-based solvents or cleaning agents.

PREVENTATIVE MAINTENANCE

Use only DeVilbiss concentrator replacement parts and accessories.

Every Drive DeVilbiss Oxygen Concentrator is tested at the factory. To assure continued trouble-free performance, the following preventative maintenance should be performed periodically by the clinician or biomedical technician. Failure to properly maintain the unit will void the warranty.

1. Check the oxygen concentration with an oxygen analyzer – every 3 years.
 - a. Calibrate the oxygen analyzer prior to checking the oxygen concentration. The analyzer should be properly calibrated using the manufacturer's recommended procedure.

NOTE– Changes in temperature, altitude or humidity may affect the analyzer's oxygen concentration reading. The analyzer should be calibrated in similar conditions to the location of the concentrator.

- b. Power the unit. Set the flow meter to 10 LPM and connect the analyzer to the unit's oxygen outlet port.
 - c. Allow the oxygen percentage to stabilize before taking a final reading, up to a maximum of 20 minutes.
 - d. Record the reading.
2. Check the audible alert and indicator lights during every service. When the power is turned ON, listen for the audible alert and check to see if the control panel LEDs are operating.
3. Inspect cabinet air filter (part #1060-608) every PM check. Replace if filter is damaged.

4. Inspect compressor intake filter (part #1060-605) at every PM check. Replace if the filter looks dirty or there is a drop in oxygen purity.
 - a. Open the filter door and replace filter as required.

5. Inspect the final bacteria filter (part #PV5LD-651) during every compressor service. Replace if dirty or there is low flow or low outlet pressure.

- a. Use the Cabinet Removal instructions found under CABINET REMOVAL in this manual to remove and attach the cabinets.
 - b. Remove the hose from each end of the filter and discard the filter.
 - c. Install the new final bacteria filter with the "IN" fitting toward the flow meter.
6. Inspect the compressor filter (part #1025D-682) during every compressor service. Replace filter if necessary.

7. Inspect the AC power cord and power switch between every patient change. Replace any damaged or defective components.

NOTE– This PM Schedule reflects:

- 4000 hour usage equal to one year
- a normal, clean operating environment.

The clinician or biomedical technician is responsible for:

- determining the condition of the concentrator operating environment.
- determining a preventative maintenance interval frequency* which takes into consideration the specific operating environment.

* Standard intervals are noted below. Service interval may be more or less frequent than stated below provided that the clinician or biomedical technician establishes and documents appropriate protocols.

PREVENTATIVE MAINTENANCE SUMMARY

Patient / Clinician

Clean and replace oxygen tubing, cannula / mask, and humidifier bottle (if used) according to manufacturer's instructions.

Clinician or Biomedical Technician

During each inspection

- Wash/Replace cabinet filter.
- Check audible alert and indicator lights.
- Inspect compressor intake filter and replace if needed

During each PM check – every 3 years for the 1060AW series

- Inspect/Replace intake bacteria filter as necessary.
- Check oxygen purity.

During compressor service

- Inspect/ Replace final bacteria filter.
- Inspect/ Replace the compressor filter

NOTE– There is no portion of the gas pathways through the concentrator that could be contaminated with body fluids under normal conditions. The device patient connection may unintentionally become contaminated with expired gases if a hose internal to the device becomes disconnected. This condition will cause no flow out of the device and/or an alarm condition.

1. Should this occur, remove the front cabinet in order to determine where the disconnection occurred.
2. Replace all components from the free end of the disconnect through the outlet port. Reference the pneumatic diagram and replacement part numbers/instructions.

MAINTENANCE

Cleaning and Disinfection When There is a Patient Change

NOTE– Recommendations for preventative maintenance at 3-year intervals are outlined in the Service and Maintenance Guidance below.

Drive DeVilbiss Healthcare recommends that at least the following procedures be carried out between uses by different patients.

NOTE– If the following described complete processing of the concentrator by an appropriately trained individual is not possible, the device should not be used by another patient.

NOTE– If preventive maintenance is due at this time, these procedures should be carried out in addition to the servicing procedures.

1. Use disinfectants safely. Always read the label and product information before use.
2. Always wear personal protective equipment when performing this procedure. Use suitable gloves and safety glasses. Cover exposed skin on arms to prevent accidental contact with bleach solution that has been applied to the concentrator.
3. Dispose of all accessories that are not suitable for reuse. This includes but may not be limited to the oxygen tubing, tubing connectors, nasal cannula and/or mask, oxygen outlet connector, and humidifier bottle. Reusing accessories marked for single use may cause risk of cross contamination or loss of performance.
4. Clean the exterior of the concentrator with a clean lint-free cloth. Heavy soil should be removed with a clean lint-free cloth dampened with water. A soft bristled brush dampened with water can be used to remove stubborn soil. Dry the concentrator using a clean lint-free cloth if water was used to remove soil.
5. Use 5.25% chlorine Sodium Hypochlorite (Clorox Regular Liquid Bleach or equivalent). Mix one (1) part Sodium Hypochlorite with four (4) parts water in an appropriate clean container. This ratio produces a one (1) part Sodium Hypochlorite to five (5) total parts solution (1:5). The total volume (amount) of solution required is determined by the number of concentrators in need of disinfection.
6. Apply the bleach solution in an even manner to the cabinet and power cord using a clean lint-free cloth. The cloth should be dampened only and not dripping of solution. Do not use a spray bottle to apply the solution. Do not saturate the device with the solution. Take care that no solution enters the vent areas on the concentrator base or door area on the back of the unit. Avoid over-saturating the cabinet seams so that no solution residue builds up in these areas. Avoid the caster wells located on the bottom of the unit.
7. Exposure time of the disinfectant solution should be 10 minutes minimum to 15 minutes maximum.
8. After the recommended exposure time, all surfaces of the concentrator should be wiped with a clean lint-free cloth dampened with drinking quality water no warmer than room temperature. Dry the unit with a dry, clean lint-free cloth. This is to remove residue that may stain or leave a film on the unit, especially after repeated disinfections.
9. Check the cord, the plug on the back of the device, the power switch, the fuse holder, and the indicator LEDs for possible damage. Replace all damaged or worn components.
10. Inspect Cabinet Air Filter and wash as needed. Replace if torn or damaged.
11. Check the oxygen concentration. If the device is within specification, the extended life intake bacteria filter does not need to be replaced between patients. If the oxygen concentration is not within specification, the provider should refer to the service manual section on Troubleshooting.
12. **OPTIONAL INSIDE CLEANING:** The concentrator must be disconnected from the power supply for this step ONLY if the unit cabinet is being opened for service: Open the concentrator and remove all dust deposits inside the

cabinet with an appropriate vacuum cleaner. Close the concentrator.

NOTE– There is no portion of the gas pathways through the concentrator that could be contaminated with body fluids under normal conditions. The device patient connection may unintentionally become contaminated with expired gases for a single fault condition i.e., a hose internal to the device becomes disconnected. This condition will cause no flow out of the device and/or an alarm condition. Should this occur, refer to the service manual for additional instructions.

Cleaning

	Recommended cleaning interval	Number of cleaning cycles *	Compatible cleaning method
Outer Cabinet	7 days	260	Water, use only a damp cloth
Filter Door Vents	7 days	260	Wipe with dry cloth, or a cloth dampened with water to remove dust.
Oxygen Outlet Connector	7 days	104	Mild dish soap (2 tbsp) (30ml) and warm water (2 cups)(400 ml)
Cabinet Air Filter	7 days inspect, clean or replace as needed	104	Mild dish soap (2 tbsp) (30ml) and warm water (2 cups) (400 ml), towel dry
Compressor Intake Filter	14 days inspect, replace as needed	---	Do not clean, replace when discolored to dark gray.

* number of cleaning cycles over the expected service life

Disinfection

NOTE– The disinfection process can only be completed by the manufacturer or by an appropriately trained technician.

	Recommended disinfection interval	Number of disinfection cycles *	Compatible disinfection method
Cabinet, power cord	Once per week when used with multiple patients.	260	1:5 chlorine Sodium Hypochlorite (5.25%)
Oxygen tubing, tubing connectors, nasal cannula/ mask, oxygen outlet connector, humidifier bottle, cabinet air filter	Do not clean, replace between patients	N/A	N/A

MAINTENANCE

SERVICE AND MAINTENANCE GUIDANCE

PULMO2 10 Liter Oxygen Concentrator Preventive Maintenance/Service Guide							
Model	Oxygen Purity Verification	Compressor Intake Filter	Internal Compressor Filter **	Final HEPA Filter **	Cabinet Air Filter	Sieve Beds **	Compressor Cup Seals **
1060AW	Every 3-years or between patients, whichever comes first or with multiple patients, every 6 months	Inspect when Recurring Reminder illuminates on control panel. Replace if filter is dirty or damaged.	Inspect in conjunction with compressor service. Replace if filter is dirty or damaged.	Inspect in conjunction with compressor service. Replace if filter is dirty or damaged.	Inspect when Recurring Reminder illuminates on control panel and wash as needed. Replace if torn or damaged.	When indicated by device performance below specification for oxygen purity, operating pressures and/or other indications of component wear	When indicated by device performance below specification for oxygen purity, operating pressures and/or other indications of component wear

** Sieve bed, compressor cup seal, internal compressor filter and final HEPA filter service should only be performed by appropriately trained biomedical technicians.

NOTE – This is a suggested maintenance and service schedule. Individual maintenance requirements may vary based upon local operating conditions, regulations, or other circumstances.

Service and maintenance should only be performed by appropriately trained biomedical technicians.

Between Patient Uses

NOTE – Scheduled maintenance should be performed in accordance with the Preventive Maintenance/Service Guide table above or Between Patient Uses.

1. Discard all oxygen tubing, cannula/mask, oxygen outlet connector, humidifier bottle, and cabinet air filter.
2. Replace cabinet air filter (when applicable) and follow the Cleaning and Disinfection Instructions.
3. Clean the concentrator cabinet and inspect/replace filters in accordance with the table above.
4. Inspect all plugs, cords, and components. Replace any damaged or worn components.
5. Check oxygen concentration with a calibrated oxygen analyzer and record the oxygen percentage. If the concentration is not within specification, refer to troubleshooting section of the IFU or the Service Manual.
6. Record the unit hours of use.
7. Verify Audible Alert and Indicator Lights at each service at startup and while operating.
8. With unit still running, unplug to test the power fail alarm. If the power fail alarm does not provide an audible alert, refer to the Service Manual or contact DeVilbiss for assistance.

SYSTEM OPERATION

The Drive DeVilbiss Oxygen Concentrator uses a pressure swing adsorption system. The air is drawn into the unit through air filters and into a double-head compressor.

A pneumatic diagram of the system is shown on page 31.

The compressed air passes through a shut-off valve before entering the rotary valve. The rotary valve is cycled at a predetermined rate, and is directed into one of two sieve beds. The sieve beds contain molecular sieve material, which is a man-made inorganic silicate. It is very porous and has the unique ability to selectively adsorb, or temporarily hold on to, nitrogen from the air as it passes through the sieve bed.

As one bed is being pressurized, the other bed is quickly depressurized. This allows the nitrogen that was adsorbed during its pressurization cycle to be exhausted from the sieve material.

The nitrogen is released through an exhaust port located on the rotary valve assembly. The nitrogen passes through a check valve before exiting through the exhaust muffler.

Also during each bed pressurization, a small amount of oxygen flows through an orifice from the pressurized bed into the depressurizing bed. This helps to clean out the nitrogen from the depressurizing bed.

The beds will continue to be alternately pressurized and depressurized as the unit operates.

Oxygen leaving the sieve beds is directed through a check valve and into the accumulator tank. A pressure regulator on the tank controls the oxygen pressure as it leaves the accumulator tank and passes through the OSD (Oxygen Sensing Device), which monitors oxygen purity and flow, before continuing to the flow meter. The flow meter allows the oxygen flow to be controlled and adjusted to the level prescribed by the patient's physician. From the flow meter, the oxygen passes through the final bacteria filter and finally the oxygen outlet port to the patient.

The Drive DeVilbiss Oxygen Concentrator operates on a timed cycle (xx sec. @ flows > xx LPM) that is controlled by the PC board. The PC board will send voltage to the valve causing it to shift and alternately pressurize the sieve beds. The PC board also activates the electronic alert system. Low flow, system abnormality, and power failure are indicated by audible and visible alerts. A high pressure condition will be indicated with a "popping" type sound produced by release of pressure from a pressure relief valve on the compressor head.

NORMAL OPERATING SEQUENCE

When the concentrator is turned ON, the following cycling sequence can be observed by attaching a pressure gauge to the accumulator tank test point.





1. The rotary valve is quickly cycled several times to relieve any residual bed pressure. This rapid cycling only happens on start-up and is clearly heard as pressure is being quickly exhausted several times. The pressure exhausts through an exhaust muffler that is connected to the valve.
2. The PC board applies a short DC voltage signal to the valve. The valve will stop for several seconds causing the right bed to pressurize first while the left bed depressurizes.
3. Voltage is again applied to the valve for a short time. The valve will stop for less than a second. During this time the sieve bed pressures are equalized. This happens too quickly to be seen on a pressure gauge.
4. A short DC voltage signal is again applied to the valve. The valve will stop for several seconds causing the left bed to pressurize while the right bed depressurizes.
5. A short DC voltage signal is again applied to the valve. The valve will stop for less than a second. During this time, the sieve bed pressures are equalized. This happens too quickly to be seen on a pressure gauge.
6. The cycle then repeats with step 2 and continues the entire time the unit is ON.

TROUBLESHOOTING

OVERVIEW OF ALARMS AND SERVICE INDICATORS

This device contains an alarm system which monitors the state of the device and alerts of abnormal operation, loss of essential performance or failures. Alarm conditions are shown on the control panel. The alarm system functions are tested at power up by lighting all visual alarm indicators and sounding the audible alarm (beep).

All alarms are Low Priority Technical Alarms.

Alert or Alarm Condition	LED Icon	Details of Alert or Alarm Condition	Visual Alert or Alarm	Audible Alarm	Action
Power-On Test	ALL	Unit was just turned On; test of all visual/audible indicators	All LEDs ON for a few seconds	Audible indicator beeps for a few seconds	Wait for unit to complete Power On Test, up to 30 seconds.
Start-up Period		The unit has recently been started and is in start-up mode, the output flow of the oxygen is temporarily < 82%	Yellow Startup Low O2 LED ON	No audible alarm during start-up period	Wait for unit to finish start-up period, up to 15 minutes
Normal Operation		Normal device operation (O2 >= 85.0%), no errors, startup period over	Green Good O2 LED ON	OFF	Use device as desired
Check Cabinet Filter Notification		Notification – Check Cabinet Filter. LED turns ON every 7 days for 12 hours then restarts the timing cycle.	Yellow Check Cabinet Filter LED ON	OFF	Check Cabinet Filter and clean (wash/dry) or replace as necessary.
Check Compressor Intake Filter Notification		Notification – Check Compressor Intake Filter. LED turns ON every 133 days for 12 hours then restarts the timing cycle.	Yellow Check Compressor Intake Filter LED ON	OFF	Check Compressor Intake Filter and replace if necessary
Low Oxygen Output Concentration		Error – O2 below 85.0% (>=60.0% < 85.0%), startup period over.	Yellow Low O2 LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Continue to use device. Refer to Troubleshooting to try to resolve problem. If problem not resolved, contact your Biomedical Technician for assistance.
Critical Low Oxygen Output Concentration		Error – O2 below 60.0% (<60.0%), startup period over.	Red Critical Low O2 LED ON. Red Service Required LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Switch to an alternate oxygen source and contact your Biomedical Technician for assistance.
Alarm Acknowledge button Off		Alarm Acknowledge button not active (either no current alarm condition or audible alarm is actively sounding).	Alarm Acknowledge button LED OFF	If no current alarm, OFF. If alarm condition is active, audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Correct the Alarm Condition if audible alert is active (beeping).
Alarm Acknowledge button On		Alarm Acknowledge button active (symbol is illuminated indicating active alarm silenced).	Alarm Acknowledge button LED ON	OFF	Correct the Alarm Condition to deactivate audible alert and return to normal operation.
Over Temperature		Warning – O2S gas over temperature or internal over temperature	Over Temperature LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Refer to Troubleshooting to try to resolve problem. If problem not resolved, contact your Technician for assistance and to arrange for servicing of the unit.
Low Flow Rate		Warning – flow rate too low (no flow or obstruction). Flow rate below low flow threshold for 30 seconds	Low Flow LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Refer to Troubleshooting to try to resolve problem. If problem not resolved, contact your Technician for assistance and to arrange for servicing of the unit.
High Flow Rate		Warning – flow rate too high. Error due to corrupted settings, OSD communication failure, fan not running, motor error, or rotary valve error. Symbol will show red.	Red Service Required LED ON	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Refer to Troubleshooting to try to resolve problem. If problem not resolved, contact your Technician for assistance and to arrange for servicing of the unit.
Device Malfunction		The device is experiencing a malfunction that requires servicing to correct. Symbol will show red.	Red Service Required LED light is illuminated	The audible alarm is beeping intermittently. Alarm can be silenced with the Alarm Acknowledge button.	Contact your service technician for assistance and to arrange for servicing of the unit.
Power Failure		Power Failure Alarm – AC power not present, microcontroller running from a backup capacitor. Symbol will show red.	Red Service Required LED is flashing	The audible alarm is beeping intermittently. Alarm Acknowledge button is not available during a power failure alarm.	Turn unit off until electric power is restored.

TROUBLESHOOTING

LED FAULT CODES - PC BOARD

Error LED	System Failure Description	Possible Causes
1	O2 Below 85%	Sieve Beds, Rotary Valve, Compressor, Leaks, Filter
2	O2 below 60%	Sieve Beds, Rotary Valve, Compressor, Leaks, Filter
3	Low O2S Flow Error	Flow Meter Set below 1lpm - Increase Flow, Leaks
4	High O2S Flow Error	Flow Meter set above 10lpm - Decrease Flow setting
5	High O2S Gas Temperature or High Enclosure Temperature	Ambient or Room temperature too high – move to a cooler location, Fan Running Slow
6	Rotary Valve Error. The unit enables the audible alert and illuminates the Service Required LED. If this error is detected, the compressor will be shut off and the system will not attempt to rotate the valve (this pre-vents damage to the stepper motor and driver)	Rotary Valve
7	Settings Error – Device settings are corrupted at startup.	PC Board
8	O2S Communications Error	PC Board
9	Fan error – the fan is not running	Fan, Wire disconnected
10	Motor Error. If this error is detected, the compressor and rotary valve will be shut off.	Motor Controller, Compressor
11	Over Pressure Error – An over pressure condition has been detected	Sieve Beds, Rotary Valve
12	Low Pressure Error	Leaks, Rotary Valve
13	Not Used	N/A

SIMPLIFIED TROUBLESHOOTING

The key to simple troubleshooting is to recognize which type of problem exists and select the most effective approach to solving the problem.

Low Oxygen Concentration

WARNING

Electric Shock Hazard. Extra care should be taken if it is necessary to operate the unit with the cabinet removed.

WARNING

Mechanical Hazard. Keep fingers, loose clothing, etc. away when working on compressor.

Observe the pressure cycle at the oxygen tank.

1. Connect a calibrated pressure gauge to the unit's oxygen tank test point or manifold. See accumulator pressure test on page 20.
2. Set the unit's flow meter to 10 LPM.
3. Power the unit and allow it to operate for a minimum of 5 minutes before observing the pressure cycle.
4. Compare the high pressures and low pressures to those expected for the current elevation and use the following chart to find the appropriate action. All consecutive high pressures should be within 2 psi of each other and all consecutive low pressures should be within 2 psi of each other.

NOTE– For normal system pressures refer to Specifications.

NOTE– Check for leaks using a certified leak detection solution such as Snoop® or equivalent (must not contain ethylene glycol). Apply leak test solution to all fittings and hose connections with unit running. If an air leak is present, the solution will bubble. All leaks should be repaired before putting the unit back in service.

WARNING

Do not apply leak test solution to any part of the rotary valve or the main PC Board assembly.

Pressure Diagnostic Chart

If Pressure Reading Is	Do This	To This
Higher than expected (pressure relief valve may be activated)	Replace	Sieve beds (2)
Lower than expected	Inspect and replace, as needed	Filters
	Inspect and correct or replace, as needed	Tubing connections
	Troubleshoot	Compressor See Component Testing, Repair and Replacement
As expected	Inspect and correct, as needed	Leaks from sieve beds to oxygen outlet port
Dropping too low during cycle change	Replace	Check valve assembly
Uneven - high pressures & low pressures are not consistent during cycle	Inspect and replace, as needed	Rotary valve wire harness
	Troubleshoot or replace	Rotary valve or PC board. See Component Testing, Repair and Replacement

TROUBLESHOOTING

TROUBLESHOOTING CHART

SYMPTOMS	POSSIBLE CAUSE	POSSIBLE REMEDY
Pulsating air noise	Intake filter not in place or defective	Check filter and replace if necessary
	Compressor intake hose disconnected	Reconnect hose
Excessive noise	Loose or defective motor mounts	Replace motor mounts
	Mounting plate loose	Reinstall or replace mounting plates
	Defective compressor	Replace compressor
	Defective cooling fan	Replace cooling fan
Fluctuating oxygen flow	Occluded humidifier	Clean or replace humidifier
	Use of improper humidifier	Use only a bubble-type humidifier designed for use with flows up to 10 liters per minute and 8.5 psi pressure
	Occluded filters	Clean or replace filters
	Occluded or defective cannula and tubing	Detach cannula from oxygen delivery tubing. If proper flow is not attained, check tubing for kinks or other obstructions. Clean or straighten as required or replace tubing if necessary
	Use of excess oxygen tubing	The unit is designed to deliver 10 lpm with a cannula on 50 feet (15 meters) of approximately 5/32" (4 mm) inside diameter tubing. Smaller diameter tubing or the addition of any other flow restriction may prevent obtaining the desired flow rate.
	Defective flow meter	Replace flow meter
	Leak in system	Check for leaks in all hoses and fittings
	Defective compressor	Replace compressor
	Defective compressor reed valve	Replace compressor reed valve
	Defective check valve	Replace check valve
Pressure regulator not adjusted properly or defective	Adjust regulator to 8.5 psi or replace regulator if it cannot be adjusted properly	
Little or no oxygen flow	Flow meter not adjusted properly	Adjust flow meter
	Hose disconnected to flow meter	Reconnect hose
	Oxygen delivery tubing is kinked or blocked	Straighten tubing or remove obstruction
	Occluded humidifier	Clean or replace humidifier
Low oxygen concentration	Leak in system	Check for leaks in all hoses and fittings
	Defective sieve bed check valve	Replace check valve
	Defective compressor reed valve	Replace compressor reed valve
	Defective compressor	Replace compressor
	Rotary valve not operating correctly	Replace valve
	Occluded filters	Clean or replace filters
	Contaminated sieve beds	Replace sieve beds
	Worn compressor cup seals	Replace compressor cup seals
Audible alarm does not sound during power failure	Unit has not been used for an extended period of time. NOTE —If the concentrator has been unused for an extended period, the unit must run for a minimum of 2 minutes before the power fail alarm will activate.	Allow unit to run for 20 minutes and retry
	Defective PC board	Replace PC board
	Defective power switch	Replace power switch
	Defective wire harness	Replace wire harness
Audible alarm does not sound when unit is turned "On"	Defective PC board	Replace PC board
Pressure relief valve activated "popping" sound	PC board connectors not properly latched	Be sure tabs are pushed completely into place
	Defective PC board	Replace PC board
	Defective rotary valve	Replace valve
	Defective rotary valve wire harness	Replace wire harness
	Defective shut-off valve	Replace shut off valve
LED(s) do not illuminate when unit is turned "On"	PC board connectors not properly latched	Be sure tabs are pushed completely into place
	Defective PC board	Replace PC board
	Defective control panel	Replace control panel

TROUBLESHOOTING

SYMPTOMS	POSSIBLE CAUSE	POSSIBLE REMEDY
No power	Line cord not properly installed or defective	Check that the power cord is properly attached to the IEC connector and is fully inserted into the receptacle. If still does not work, replace power cord.
	No power at receptacle	Check building circuit breaker or fuse, or have house wiring checked by qualified electrician. Circuit may be fully loaded with other appliances and another receptacle may be required.
	Defective power switch	Replace power switch
	Defective fuses	Replace fuses
	Defective controller	Replace controller
Compressor is not running, but fan is	Main wiring harness disconnected/defective	Reconnect/replace wiring harness
	Loose compressor wire	Tighten or attach wire
	Defective compressor	Replace compressor
Unit warm to the touch and cannot be restarted for several minutes	Compressor overheated due to: 1. Occluded filters 2. Restricted input or output air passage 3. Low or high line voltage	1. Clean or replace filters 2. Remove obstruction 3. Check line voltage; use alternate circuit independent of other appliances
	Defective cooling fan	Replace cooling fan
	Defective compressor	Replace compressor
Fluctuating or no flow	System pressure below 8.5 psi (59 kPa) due to:	
	1. Leak in system	1. Check for leaks in all hoses and fittings
	2. Defective compressor	2. Replace compressor

NOTE-If unit does not run or perform to specification, verify the proper AC input voltage and frequency are being used. The AC input voltage should be between 90 and 264 volts

PROPER REPAIR PROCEDURES

WARNING

When servicing the Drive DeVilbiss Oxygen Concentrator, be absolutely certain that the correct tools are used and that the parts are free of oil and grease or any material not compatible with oxygen. Teflon® tape is recommended and must be applied to the male threads omitting the first thread to eliminate the possibility of tape particles entering the oxygen system.

Electric shock hazard. Do not remove cabinet. The cabinet should only be removed by a qualified Drive DeVilbiss homecare provider.

Disconnect the power cord from the wall outlet before attempting repairs on the unit. Extra care should be taken if it is necessary to operate the unit with the cabinet removed.

NOTE- Be sure to read all of the steps involved before beginning any of the procedures in this manual.

NOTE- After repairing or replacing a component check the oxygen concentration and test for leaks.

The Drive DeVilbiss Oxygen Concentrator is designed for ease of service. To aid service personnel a Service Kit (part #444-501) is available which contains the necessary gauges, tools, and testing instruments to properly service the oxygen concentrator. See list below.

In addition, you will also need an oxygen analyzer (part #R217P62) to periodically check oxygen concentration levels and leak test solution.

The following parts are included in the Service Kit:

1	Slotted bit
1	#1 Phillips bit
1	#2 Phillips bit
1	Crescent wrench
1	8" Duckbill pliers
1	Voltmeter
2	Pressure/Vacuum gauge
1	Tool box
2	Test Fittings
1	Torx screwdriver w/bits
1	Channel Lock Pliers
1	1/4" Ratchet wrench
1	10mm Socket 1/4" Drive
1	1/4" Drive extension
1	Plastic storage case

CABINET REMOVAL

The majority of all the servicing and repairs can be done without removing the front bib. However, to gain access to the components behind the bib, it may be loosened or removed.

To remove back cabinet:

1. Ensure the unit is unplugged from the wall outlet.

COMPONENT TESTING, REPAIR AND REPLACEMENT

2. Remove the six screws that hold the back cabinet to the internal structure and the bib - 2 behind the handle, 2 at the bottom and 1 on each side.
3. Remove the back cabinet by sliding it toward the rear until clear.

To remove the front cabinet:

4. Remove the front cabinet by pushing the top shoulders toward the back of the unit, then outward away from behind the bib. Tilt the top of the front cabinet forward until it can be pulled out of the base of the unit.

To loosen the bib:


1. Remove the two screws (located directly above the hour meter) that hold the bib to the unit's internal structure. This will allow access to the components behind the bib.

To remove the bib completely:

1. Remove the two screws as above.
2. Disconnect the ribbon cable at the connector.
3. Unplug the power switch connector from the main wire harness.
4. Tilt the top of the bib forward to release it from the slot in the body of the concentrator.
5. Remove the hose connected to the bottom of the flow meter.

To reassemble bib:

1. Reconnect the wires and hose.
2. Insert the bib tab into the slot above the Rotary valve, and push until it snaps into place.
3. Secure bib with two screws.

Accumulator "T" Fitting

To remove the accumulator tank:

1. Remove the blue tubing from the accumulator tank "T" fitting and the pressure regulator.
2. Loosen the ladder clamp that secures the accumulator tank to the unit.
3. Slide the tank out and away from the unit to remove.


ACCUMULATOR TANK

The oxygen accumulator tank holds the concentrated oxygen and releases it to the patient at a specified liter flow.

ACCUMULATOR PRESSURE TEST

To check accumulator pressures:

1. Make sure the unit is "Off."
2. Use the Cabinet Removal instructions listed previously to open the unit for testing.
3. Use the pressure gauge (part #PVO2D-601) and pressure test assembly (part #303DZ-637) included in the Service Kit.

4. Remove the tubing cap from the accumulator tank fitting attached to the tank, and attach the 1/16" (1.6 mm) diameter tubing from the gauge to the fitting just vacated above.

See picture below showing pressures being checked at accumulator tank "T" fitting.

COMPONENT TESTING, REPAIR AND REPLACEMENT

- Turn the unit "On" with the flow rate set to maximum recommended flow, which is 10 lpm. Allow the unit to run for 5 minutes before observing the pressures. During each timed cycle, the average pressure in the oxygen accumulator will rise and fall. The high pressures should be consistent and the low pressures should be consistent. The pressure swing will be approximately 4-5 psi.

NOTE– Expected normal pressures observed depend on altitude and flow rate. See the Typical Peak Accumulator Tank Pressure Range chart below.

- Increases in altitude and flow rate will slightly decrease accumulator pressures.
- Lower altitudes and flow rates will slightly increase accumulator pressures.

NOTE– A defective check valve in the purge harness may cause a rapid drop in accumulator pressure below the minimum value.

TYPICAL PEAK ACCUMULATOR TANK PRESSURE RANGE @ 10LPM		
Altitude	Psi	kPa
0 to 457 m 0 to 1500 ft.	25-36	172-248
457 to 914 m 1500 to 3000 ft.	21-33	145-228
914 to 1524 m 3000 to 5000 ft.	21-30	145-207

- Refer to the Type 1 – Purity Issues, found under Simplified Troubleshooting, to determine the appropriate action to take in resolving abnormal pressure cycles.

NOTE– A defective compressor will be indicated by slowly rising pressure. Pressure may only reach a certain level and then stop.

Low oxygen concentration levels and accumulator pressures higher than normal may indicate defective sieve beds. Severely contaminated beds may also cause the pressure relief valve on the compressor to open.

NOTE– A malfunctioning rotary valve may also cause high accumulator tank pressure and activation of the pressure relief valve. In this case it should be determined whether the problem is with the sieve beds, valve, or both.

CHECK VALVES


Sieve Bed Check Valves

The sieve bed check valves, located in the check valve tubing assembly, between the outlet of each sieve bed and the accumulator tank, allow oxygen to pass from the sieve beds to the accumulator tank when the bed pressure is greater than the accumulator tank pressure. These valves also prevent reverse flow of oxygen from the accumulator tank to the sieve beds. A defective sieve bed check valve will result in lower oxygen concentrations and accumulator pressures.

A fixed orifice is located in the orifice tube assembly and is connected between the two sieve beds. This orifice directs a small amount of pressurized oxygen into the discharging sieve bed to help with the nitrogen exhaust process.

Use the Accumulator Pressure Test listed previously to troubleshoot sieve bed check valves.

- To replace defective sieve bed check valves, remove the tubing on from the top of each sieve bed and also the tubing that connects to the T fitting on the accumulator tank, then install new check valve tubing assembly making sure the arrow is pointing away from the sieve bed.



Exhaust Check Valve

The exhaust check valve, located between the rotary valve and the exhaust muffler, prevents room air from entering the sieve beds, thus preventing sieve material contamination and lengthening component life.

A defective exhaust check valve may cause high pressure in the unit.

COMPRESSOR

The Drive DeVilbiss Oxygen Concentrator uses a double-head, oil-free compressor. The compressor is powered by a brushless DC motor. It is secured to the base with a mounting plate and four motor mount springs.

Sleeve w/Heat Sink

Valve Plate

A compressor that is worn or defective may:




- cause pressure to rise slowly.
- cause excessive noise and/or vibration.
- cause lower oxygen concentrations.

A worn or defective compressor can be caused by a defective internal component such as:

- reed valve
- o-ring
- gasket
- cup seal

These components are included in the Compressor Rebuild Kit (1025K-643).

COMPONENT TESTING, REPAIR AND REPLACEMENT

NOTE— A built-in thermal cutoff switch will shut the compressor off if it becomes overheated. This protects the compressor from damage caused by heat build-up. (Some models have an auxiliary thermostat mounted within the compressor compartment.) Should this condition occur, the compressor will require several minutes for the thermo-protective device to reset.

NOTE— A pressure relief (PR) valve is located on the pressure head to prevent high pressure build up in the system should a component malfunction occur.

To test the compressor for proper output:

NOTE— If the compressor is not providing a high enough output the patient alert system may be activated.

1. Use the Cabinet Removal instructions listed previously to open the unit for testing.
2. Use the Accumulator Pressure Test listed previously to observe the high pressures, low pressures and the pressure drop.
3. Refer to the Low Oxygen Concentration, found under Simplified Troubleshooting, to determine the appropriate action to take in resolving abnormal pressure cycles.

NOTE— A compressor, which slowly builds pressure that remains below 25 psi, indicates worn cup seals and/ or reed valves.

If these conditions are observed then:

- The unit filter(s) may be occluded—check the air filter, compressor filter, and compressor intake filter for occlusions.
- There may be a severe leak in the system—check for air leaks using a leak detection solution such as Snoop[®] or equivalent (must not contain ethylene glycol).

CAUTION

Do not apply leak test solution to any part of the rotary valve or the main PC Board assembly.

- The compressor reed valves, cup seal, or the compressor itself may be defective.

If the filters are not occluded and no leaks are found, the compressor must then be removed and repaired or replaced.

COMPONENT TESTING, REPAIR AND REPLACEMENT

To remove the compressor:

1. Make sure the unit is unplugged from the wall outlet.
2. Use the Cabinet Removal instructions listed previously to open the unit.


CAUTION

If the unit has been running recently, the compressor may be hot.


3. Disconnect the connector on the compressor wire harness from the controller and remove harness from compressor box.
4. Disconnect the intake hose from the bottom of the grommet.
5. Remove the compressor filter by first, loosening the ladder clamp on the heat exchange tube and push down to disconnect. This can be done using a flat head screw driver or a 1/4" nut driver. Next, loosen the ladder clamp on the exhaust fitting and remove filter. You may need to rotate the exhaust fitting down a bit to remove the hose from the fitting.
6. Remove 2 screws from front of compressor retaining plate.
7. Remove 2 hex nuts from the other side of the unit located in recesses on either side of the rotary valve.
8. Carefully lift the compressor and mounting plate assembly out of the compressor box.

If the compressor will be replaced:

9. Disconnect the intake hose from both intake fittings on compressor.

10. Turn the compressor upside down and place on work surface.
11. Remove both retaining plates on each side of compressor by twisting plates away from the unit until free of mounting plate.


12. Remove springs and black screws/bumpers holding the mounting plate to the compressor.

To inspect and/or replace internal components:

1. Place compressor upright and remove the eight screws that hold the compressor heads in place. When removing the heads, be sure to keep

each head and its components with the correct compressor side.


2. Check for proper placement of or damage to the gaskets on the bottom of the compressor heads. Replace if damaged.

3. Remove reed valve plates. A reed valve is located on each side of the valve plate.
4. The compressor reed valves should be flush with the valve plate. If the valve is broken or not flush with the valve plate, or foreign matter is detected inside the head, clean or replace the compressor reed valves.

To replace the compressor reed valves:

- a. Remove the screw holding the compressor reed valves in position on the valve plate and discard the used reed valves.
 - b. Position the new reed valves so that they are centered and completely cover the holes in the valve plate.
 - c. Place the metal retainers on the reed valves and secure with the reed valve screw.
5. Check for proper placement of or damage to the rubber o-ring on the bottom of the valve plate. Replace if damaged. Refer to the compressor rebuild kit pictures on page 22.
 6. Remove piston sleeves by pulling upward and inspect cup seal on pistons. Replace if badly worn or damaged.

To replace cup seal:


- a. Remove both rod screws from top of piston. Note the position of the screws since they are different lengths.
 - b. Remove the cup retainer plate.
 - c. Discard defective cup seal.
 - d. Place new cup seal into position.
 - e. Replace cup retainer plate.
 - f. Secure with screws.
7. Clean inside surface of sleeves before reinstalling. Position sleeve at 45 degree angle over the piston. Carefully push it down as you rotate it slightly around the top of the piston until it is in place.

COMPONENT TESTING, REPAIR AND REPLACEMENT

- Place valve plates on the compressor so that heads of reed valve screws are aligned with the indentation in top of pistons.
- Install the compressor heads so that the holes in the heads are aligned with the holes in the compressor housing.
- Secure compressor heads with the screws.

To replace the compressor:

- Position mounting plates with hooks facing intake side of compressor and secure with four bumpers and screws.

- Install the compressor mounting springs over the bumpers and carefully reattach the retaining plates, making sure that the springs stay between the bumpers on both plates.

- Turn the compressor upright and attach the intake hose to both intake fittings and secure with ladder clamps.
- Carefully place the compressor in the compressor box so that the studs on the back of the retaining plates line up with the screw holes on the unit.

NOTE– Ensure the compressor mounting and retaining plates are positioned properly to prevent excessive noise and vibration.

- Reinstall the 2 hex nuts onto the studs on the back of the retaining plates through the recesses on either side of the rotary valve.
- Reinstall the 2 screws on the front of the compressor retaining plate.
- Reconnect intake hose to the bottom of the compressor intake filter grommet.
- Route compressor wire harness under the foam, then under the retaining plate on the left. Then run the harness through the hole located on the side of the compressor box near the bottom. Run the harness up the side of the compressor box, pushing the harness behind the two tabs on the side. Secure harness and push the connector securely into the MOTOR OUT connector on the controller.
- Install compressor filter by attaching tubing and ladder clamps to compressor exhaust fitting and heat exchange tube. Ensure heat exchange tube is not touching mounting plate or fan.


CONTROLLER (POWER SUPPLY/CONTROL MODULE)

The controller uses a switched-mode power supply to convert the incoming AC voltage to a DC output.

The controller also uses a high voltage brushless DC motor driver to control the compressor. The motor driver can read and set the compressor speed as well as read the voltage, current, temperature and error flags.

To test the power supply:

- Make sure the controller is getting power by checking the AC voltage at the connector of the main wire harness that connects to the controller.
- If the unit is getting appropriate incoming AC voltage, check the DC output voltage coming from the DC OUT connector on the controller.

- NOTE** - Use caution when connecting a multimeter. Do not short terminals.
- If the controller is not putting out 12 volts DC, replace the controller.

To replace the controller:

- Disconnect the unit from AC power.
- Follow instructions for removing the cabinet and bib.
- Remove all wires connected to the controller.
- Remove the 2 screws at the top and 2 screws at the bottom of the controller that are securing it to the mounting bracket.
- Replace with new controller and secure to the mounting bracket using the screws previously removed.
- Reconnect wires.

COOLING FAN

The cooling fan provides constant air flow to cool the compressor. The cooling fan is located in the bottom of the unit, below the compressor.

The controller will self-adjust the fan speed based on internal temperature. If the temperature reaches 75°C +/-5; the thermal cut off switch, located in the compressor box, will cut power to the unit to prevent damage. Once the unit cools down to 60°C +/-10 the thermal cut off switch will reset, allowing the unit to be powered on.

To replace the cooling fan:

- Make sure the unit is unplugged from AC power.
- Use the Cabinet Removal instructions listed previously to open the unit.
- Use the To Remove the Compressor instructions listed under Compressor.
- Disconnect the cooling fan connector from the PC board.
- Note the position of the fan and wires before removing the four retaining screws that secure the fan to the base of the unit.
- Remove the defective fan and secure the replacement fan in position with the four retaining screws.

NOTE - When installing the fan, be sure the air flow directional arrow on the side of the fan is directed away from the compressor and fan guard is reinstalled properly.

- Reconnect the electrical connector.
- Reinstall the compressor.

FLOW METER

The DeVilbiss 1060AW series oxygen concentrator flow meter (1060-607) has an operating flow rate of 1-10 LPM. Flows below 1 LPM may cause the low flow alarm to activate.

Locking Screw

The flow meter is pressure compensated and has an accuracy level of $\pm 5\%$ of full scale at all liter flows except at 10 LPM. Accuracy at 10LPM is $+0\%$ / -5% . The flow meter can be locked using a 1/16th inch Allen wrench and tightening the locking screw behind the flow meter knob.

To check for leaks in the flow meter tubing:

1. Check for leaks using a certified leak detection solution such as Snoop® or equivalent (must not contain ethylene glycol).
2. Apply leak test solution to all fittings and hose connections with the unit running.

CAUTION

Do not apply leak test solution to any part of the rotary valve or the main PC Board assembly.

3. If an air leak is present, the solution will bubble. All leaks should be repaired before putting the concentrator back in service.

WARNING

Electric Shock Hazard. Use caution when leak testing near electrical connections.

To replace the flow meter:

1. Make sure the unit is unplugged from the wall outlet.
2. Remove the front and back cabinets.
3. From behind the bib, remove the 2 hoses from the flow meter.
4. Unscrew both the top and bottom flow meter fittings and remove the flow meter by pushing it out through the bib.
5. Install new flow meter and reconnect hoses.

FUSED IEC CONNECTOR

The IEC connector for the power cord contains two AC fuses. These fuses provide overcurrent protection to the unit. In the event of an overcurrent condition, the fuses will open and cut off power to the unit. To check the fuses, first make sure that the unit is disconnected from AC power and remove the power cord from the IEC connector. Slide open the fuse holder by pulling it forward and push out the fuses.

To test the fuse using continuity, set the digital multimeter to measure continuity. Place each probe on opposite ends of the fuse (polarity does not matter). Listen for the multimeter to beep continuously as you hold the probes against the fuse. If you don't hear any noise coming from the meter, then the fuse is open and should be replaced.

To test the fuse using resistance, set the digital multimeter to measure resistance (Ohms). Touch the probes together to get an initial reading. Then put the probes on either side of the fuse and check if the reading is similar. If it is, then the fuse works properly. If you get no reading or "OL", then the fuse is open and should be replaced.

HOUR METER

To replace the hour meter:

1. Make sure the unit is unplugged from the wall outlet.
2. Remove the front and back cabinets and loosen the bib.
3. Disconnect the hour meter connector from the PC board.
4. Remove the meter by carefully inserting a small flat screwdriver under the outer edge of meter and prying upward.
5. Install a new hour meter by applying downward pressure until it snaps into position.
6. Connect the hour meter to the PC board.

LEAK TESTING

1. Check for leaks using a certified leak detection solution such as Snoop® or equivalent (must not contain ethylene glycol).
2. Apply leak test solution to all fittings and hose connections with the unit running.

CAUTION

Do not apply leak test solution to any part of the rotary valve or the main PC Board assembly.

3. If an air leak is present, the solution will bubble. All leaks should be repaired before putting the concentrator back in service.

WARNING


Electric Shock Hazard. Use caution when leak testing near electrical connections.

NOTE– Additional leak testing can be done by occluding the outlet port with a finger or thumb while the unit is running. When the port is occluded the flow meter ball should drop to zero. If flow still registers on the flow meter, this indicates a leak. Check for leaks around the flow meter, final bacteria filter and outlet port.

MOLECULAR SIEVE BEDS

The build and release of pressure in the sieve beds indicates the health of the sieve material and the operation of contributing components. Determine 'good' or expected pressures for your altitude by testing bed pressures on multiple oxygen concentrators that are producing at least 93% purity.

Check sieve bed pressures at the accumulator tank using the accumulator tank "T" fitting. See picture below.

Accumulator "T" Fitting

NOTE– The pressure will rise and fall as it cycles through the rotary valve. The pressure swing will be approximately 4-5 psi.

Also consider the acceptable pressure range for various altitudes as shown below.

TYPICAL PEAK ACCUMULATOR TANK PRESSURE RANGE @ 10LPM

Altitude	Psi	kPa
0 to 457 m 0 to 1500 ft.	25-36	172-248
457 to 914 m 1500 to 3000 ft.	21-33	145-228
914 to 1524 m 3000 to 5000 ft.	21-30	145.207

To test sieve bed pressures:

- Remove the plugged piece of 1/16" tubing from the accumulator tank "T" fitting.
- Connect the pressure gauge to the test point and observe the cycling of pressures. Refer to the section on Normal Operating Sequence and also Accumulator Pressure Test.
 - If it is determined that the valve is not cycling the pressure, refer to the section on Rotary Valve testing.
 - If the unit is cycling properly; allow the oxygen percentage and accumulator tank pressure to stabilize, up to a maximum of 20 minutes.
- Observe the high and low pressures and check oxygen concentration.
 - If the pressures are within the expected range and the oxygen concentration is within specification, the sieve material is effective.
 - If pressures are not within the expected range and/or the oxygen concentration is not within specification, refer to the Simplified Troubleshooting Table in this manual to determine the proper corrective action.

NOTE– If the molecular sieve material is found to be no longer effective, first search for the source of a malfunction in the system; then, for a cause for contamination (such as leaks) and take corrective action.

To replace the molecular sieve beds:

NOTE– Make sure that the sealing caps remain on the new sieve beds until just prior to connecting hoses and tubing.

- Make sure any contamination problem has been corrected before replacing.
- Make sure the unit is unplugged from the wall outlet.
- Loosen the ladder clamps that secure the sieve beds to the internal structure of the unit.
- Remove the tubing from the fittings at the top of each sieve bed.
- Remove the hose clamps and hose from the bottom of each bed.
- Lift each sieve bed upward until they are free of the ladder clamps.
- Install new sieve beds in reverse order. Make sure to correctly position and reattach hose and hose clamps to the bottom of the beds before tightening the ladder clamps.
- Leak test all connections with a certified leak detection solution such as Snoop® or equivalent (must not contain ethylene glycol). Apply leak test solution to all fittings and hose connections with unit running. If an air leak is present, the solution will bubble. All leaks should be repaired before putting the unit back in service.

CAUTION


Do not apply leak test solution to any part of the rotary valve or the main PC Board assembly.

POWER SWITCH

To replace the power switch:

- Make sure the unit is unplugged from the wall outlet.
- Remove front and back cabinets and loosen the bib.
- Disconnect the power switch connector from the main wire harness and disconnect the red wires from the PC board.
- While squeezing the locking tabs on the top and bottom of the switch, push the switch out of the front of the unit.
- Install the new switch in the correct orientation making sure that it locks into position.
- Reconnect the power switch connector to the main wire harness and connect the red wires to the PC board.


Switch Detail

*Either brown wire can be connected to either left hand terminal of the switch.
**Either red wire can be connected to either right hand terminal of the switch.

PRESSURE REGULATOR

The pressure regulator, located at the top of the accumulator tank, is connected to the OSD, and keeps the flow of oxygen to the patient at a steady 8.5 +/- 0.5 psi (58.6 +/- 3.5 KPa)

Regulator for Outlet Pressure to Patient

It is not necessary to test the pressure regulator unless there is a problem with flow rate accuracy. The following test is only needed when troubleshooting flow rate issues.

To test the pressure regulator:

- Turn the unit "On."
- Set the flow meter at 10 lpm
- Attach a pressure gauge (part #PVO2D-601) to the oxygen outlet to obtain a reference pressure. Use this reference pressure to determine if further testing is needed.
- If the reference pressure varies from the expected pressure by more than ± 1 psi or ± 7 kPa, connect a 'T' fitting directly between the pressure regulator fitting and the pressure regulator tubing and attach the pressure gauge to the 3rd leg of the 'T' fitting. If the pressure reading is not within 8.5 ± 0.5 psig (58.6 ± 3.5 kPa), adjustment to the pressure regulator is required.

COMPONENT TESTING, REPAIR AND REPLACEMENT

To adjust the pressure regulator:

1. If necessary, the pressure regulator can be adjusted by turning the allen screw on top of the regulator until the pressure is within specification. Turn clockwise to increase the pressure, and counterclockwise to decrease the pressure.

NOTE– Before adjusting the pressure regulator, make sure no leaks exist by using a certified leak detection solution such as Snoop® or equivalent (must not contain ethylene glycol).

CAUTION

Do not apply leak test solution to any part of the rotary valve or the main PC Board assembly.

A malfunction in the pressure regulator will cause either a loss or fluctuation in the oxygen flow which will be seen on the flow meter or a decrease in oxygen concentration.


To replace the pressure regulator:

1. Make sure the unit is unplugged from the wall outlet.
2. Use the Cabinet Removal instructions listed previously to open the unit.
3. Remove the tubing from the pressure regulator.
4. Loosen the ladder clamp on the accumulator tank and rotate it forward toward the front of the unit, just enough to allow the regulator to be removed.
5. Unscrew the regulator from the accumulator tank.
6. Install a new regulator on the accumulator tank and verify correct pressure using a calibrated pressure gauge. Adjust if necessary.
7. Rotate the accumulator tank back to its original position and tighten ladder clamp.
8. Reconnect tubing to the pressure regulator.

PRINTED CIRCUIT BOARD

The printed circuit (PC) board is responsible for monitoring and controlling the Drive DeVilbiss Oxygen Concentrator; it contains a series of LEDs to help diagnose problems. See the LED Fault Code chart in Troubleshooting.

The PC board has preset alerts for low flow and power failure. Should any of the alert values be exceeded, the patient alert system will activate.

NOTE– If the concentrator has been unused for an extended period, the unit must run for 2 minutes before the power fail alarm will be enabled. This alarm is powered by a capacitor on the PC board. The capacitor will be fully charged after 10 minutes.


CAUTION

Do not apply any force or flex to the PC Board when connecting or disconnecting electronic or pneumatic components. Damage to the electronic assembly is possible.

To remove and replace the PC board:

1. Make sure the unit is unplugged from the wall outlet.
2. Use the Cabinet Removal instructions listed previously to open the unit.
3. Follow instructions for removing the accumulator tank.
4. Disconnect all wires and electrical connectors.
5. Remove the 1/8" (3.2mm) tubing from both fittings on the oxygen sensor.
6. Disconnect the coupler on the 1/8" (3.2mm) tubing connected to the pressure sensor from the tubing that is going to the accumulator tank.
7. Remove the 3 screws that secure the board to the mounting plate and remove the PC board.
8. Install the new PC board and secure it using the 3 screws.
9. Reconnect all wires, electrical connectors and tubing.
10. Reinstall the accumulator tank.

PC Board Connector Diagram

Rotary Valve

The rotary valve alternately distributes compressed air, supplied by the compressor, to the sieve beds. While one bed is being pressurized the other bed is being exhausted through the valve exhaust port.

The valve contains two revolving discs powered by a stepper motor to cycle the pressure between the beds. DC voltage is supplied by the PC board to the motor windings causing the internal discs to turn and direct pressure to the proper sieve bed.

If the rotary valve is not shifting properly, one bed may pressurize continuously. The pressure can build up until the pressure relief valve on the compressor releases the excess with a sputtering sound.

There are several reasons why the rotary valve could malfunction, so the cause of failure must be determined before corrective action can be taken.

To test the rotary valve:

1. Use the Cabinet Removal instructions listed previously to open the unit for testing.

2. Connect a pressure gauge to the test point on the accumulator tank to observe unit cycling and the high and low bed pressures. Refer to the section on Normal Operating Sequence and also Accumulator Pressure Test.

- a. If it is determined that the valve is shifting properly and the high and low pressures are consistent, the valve is operating correctly.
- b. If it is determined that the valve did not shift properly or the high pressures or low pressures are uneven, check the valve wire harness for any loose or broken connections. If harness connections are good, then the problem is caused by either the rotary valve or PC board.

To replace the rotary valve:

1. Make sure the unit is unplugged from AC power.
2. Follow the cabinet removal instructions previously listed.
3. Unplug the rotary valve wire harness from the stepper motor on the valve.
4. Loosen the ladder clamps and remove the hoses that go to the sieve beds from the rotary valve.
5. Loosen the ladder clamp on the back of the valve and remove the exhaust hose from the valve.
6. Using a crescent wrench, loosen the fitting at the bottom of the rotary valve and unscrew it from the shut-off valve.
7. Screw the new rotary valve into the shut-off valve and tighten securely.
8. Reconnect hoses and wire harness.

Shut-Off Valve

The shut-off valve is designed to close when the unit is powered off. This prevents room air from entering the sieve beds, which will decrease sieve material contamination and help lengthen the component life.

To test the shut-off valve:

1. Follow the directions for removing the rotary valve.
 2. Power on the unit and verify that air is coming out of the shut-off valve.
- i** NOTE- The unit will not run for long before detecting that the rotary valve is disconnected.
3. If no air comes out of the shut-off valve, it is defective and must be replaced.

To replace the shut-off valve:

1. Make sure the unit is unplugged from AC power.
2. Follow directions for removing the rotary valve.
3. Unplug the shut-off valve wire harness from the PC board.
4. Loosen the ladder clamp on the hose connecting to the heat exchange tube and remove hose from valve fitting.
5. Remove the defective valve and unscrew the fitting using a crescent wrench.
6. Install the fitting that was just removed on to the new shut-off valve.
7. Reconnect hose and wire harness.
8. Reinstall the rotary valve.


FIGURES, DIAGRAMS AND PARTS LIST

1060AW Unit



Front View

Visible Alerts



Rear View


Filter door with cabinet gross particle air filter (and foam)


Front Inside


Intake Filter behind Filter Door

Accumulator Tank (close-up)


Accumulator Tank (removed)


FIGURES, DIAGRAMS AND PARTS LIST

1060AW Unit



Rear Inside

Behind Front Cover

Inside Compressor Box

	PART DESCRIPTION	PART NUMBER
1	Accumulator Tank	1060-610
2	Accumulator Tank T Fitting	444-582
Cabinet Parts:		
3	Base Assembly - includes IEC connector, foam, 4 fan isolation mounts, cable ties	1060-603
4	Bib	1060-613
5	Compressor Box	1060-614
6	Filter Door	1060-639
7	Front Cover	1060-611
8	Rear Cover	1060-612
9	Caster, non-locking (2 ea), Caster, locking (2 ea)	1060-650
10	Control Panel (1 ea)	1060-631
11	Cooling Fan (1 ea)	1060-634
12	Compressor	1060-625
13	Compressor Mounting Plate & Retaining Bracket - set	1060-632

14	Compressor Mounting Springs	1060-609
15	Compressor Service Kits (2 ea)	1060-643
16	Controller (Power supply / control module)	1060-621
17	Exhaust Check Valve	1060-628
18	Exhaust Muffler	1025D-705
19	Elbow, Muffler (2 ea)	1060-656
20	Elbow, Intake fitting (1 ea)	1060-657
	FILTERS:	
21	Cabinet Filter (1 ea)	1060-626
22	Compressor Filter	1025D-682
23	Filter Package - cabinet filter, intake filter, final bacteria filter, compressor filter (1 ea)	1060-627
24	Final Bacteria Filter	PV5LD-651
25	Compressor Intake Filter (2 ea)	1060-605
26	Flow Meter 10L	1060-607
27	Heat Exchange tube	1060-655
28	Hour Meter	1060-660
	HOSES:	
29	Blue Pressure Hose (34 in.)	1060-671
30	Braided Silicone Pressure Hose (16.5 in.)	1060-670
31	Compressor Intake hose	1025D-657
32	Check Valve Tubing Assembly	1060-629
33	Formed Braided Silicone Pressure Hose (set of 2)	1060-654
34	Orifice Tube assembly	1060-653
35	IEC Connector -Fused with AC Harness 5 x 20 mm 10A slo-blo	1060-624
36	Intake Filter Grommet (1 ea)	1060-604
37	Ladder clamps – sieve beds, accumulator tank	1060-642
38	Mounting Bracket for Controller / PCBA Board	1060-623
39	Oxygen Outlet Port Assembly (1 ea)	1060-602
40	PC Board with tubing	1060-622
41	Power Cord, Locking – European plug (2 ea)	1060-615
42	Power Cord, Locking – UK plug (2 ea)	1060-616
43	Power Cord Strap (10 ea)	MC29D-657
44	Power Switch (1 ea)	1060-608
45	Pressure Regulator 8.5 PSI	MC29D-612
46	PR valve (45 +/- 5 PSI)	515ADZ-614


47	Replacement fuses (20 ea) 5 x 20 mm 10A slo-blo	1060-617
48	Rotary Valve assembly	1060-640
49	Shutoff valve	1060-641
50	Screws - all screws used in the assembly	1060-644
51	Sieve Bed Set	1060-619
52	Thermostat	1060-630
	WIRE HARNESS:	
53	Main Wire Harness	1060-620
54	Rotary Valve Harness	1060-672
55	Complete Repair Parts Kit EU AC cords	1060-680EU
56	Complete Repair Parts Kit US AC cords	1060-680US
57	Complete Repair Parts Kit UK AC cords	1060-680UK
58	Power Cord, Locking- US plug (2 ea)	1060-618
59	External Filter Set; Cabinet Filter (2 ea), Intake Filter (2 ea)	1060-675

	ACCESSORIES	
	High Flow Bubble Humidifier, 6-15LPM (1 ea)	7900-0
	Oxygen Outlet Connector (1 ea)	CN100
	Packaging kit	1060-606
	Nasal Cannula Pediatric	1602-7
	Nasal Cannula Neonatal	1611-7
	High Flow Nasal Cannula	1600HF-7

	TOOLS	
	Service Kit	444-501
	Pressure Gauge	PVO2D-601
	Pressure Test Assembly	303DZ-637
	Oxygen Analyzer (Max O2)	R217P62


FIGURES, DIAGRAMS AND PARTS LIST

Pneumatic Diagram

FIGURES, DIAGRAMS AND PARTS LIST

Wiring Diagram

ORDERING INFORMATION AND PARTS RETURN

ORDERING INFORMATION

When ordering components, instruction guides, or service manuals the following must be provided:

- Unit Catalog Number
- Unit Serial Number
- Part Number
- Quantity Required

Drive DeVilbiss Concentrator 1060AW 10 Liter Series Instruction Guide:

- SE-1060AW 1060AW English, Spanish, French, Amharic, Hindi

Drive DeVilbiss Concentrator 1060AW 10 Liter Series Service Manual:

- LT-2355 1060AW English
- LT-2355-ES 1060AW Spanish
- LT-2355-FR 1060AW French
- LT-2355-AM 1060AW Amharic
- LT-2355-HI 1060AW Hindi

Orders may be placed by calling:

- Customer Service: 800-338-1988
- International Department: 814-443-4881
- Europe: +49 (0) 621-178-98-0

RETURNS

Before returning units to the factory, call the Drive DeVilbiss Healthcare Customer Service Department (800-338-1988) or (814-443-4881) to obtain a return authorization number. Include in the package a note indicating the return authorization number along with your company name, address, phone number, and account number. The return authorization number should also be written on the outside of the package.

To expedite your order for non-warranty parts, the following information should be given to the representative:

- Catalog number
- Serial number
- Hour meter reading for each concentrator
- Account number
- Company name and address
- Description of problem

WARRANTY

WARRANTY

Drive DeVilbiss Healthcare warrants the Drive DeVilbiss 1060AW 10 Liter Oxygen Concentrator under the conditions and limitations stated below. Drive DeVilbiss warrants this equipment to be free from defects in workmanship and materials for three (3) years from date of factory shipment to the original purchaser, (typically the healthcare provider) unless contractually specified otherwise. This warranty is limited to the Buyer of new equipment purchased directly from Drive DeVilbiss, or one of its Providers, Distributors, or Agents. Drive DeVilbiss' obligation under this warranty is limited to product repair (parts and labor) at its factory or at an Authorized Service Center. Routine maintenance items, such as filters, are not covered under this warranty, nor does it cover normal wear and tear.

Warranty Claims Submissions

The original purchaser must submit any warranty claim to Drive DeVilbiss or to an Authorized Service Center. Upon verification of the warranty status, instructions will be issued. For all returns, the original purchaser must (1) properly package the unit in a DeVilbiss approved shipping container, (2) properly identify the claim with the Return Authorization Number, and (3) send the shipment freight prepaid. Service under this warranty must be performed by Drive DeVilbiss and/or an Authorized Service Center.

i NOTE– *This warranty does not obligate Drive DeVilbiss to provide a loaner unit during the time that an oxygen concentrator is undergoing repair.*

i NOTE– *Replacement components are warranted for the unexpired portion of the original Limited Warranty.*

This warranty shall be voided, and Drive DeVilbiss shall be relieved of any obligation or liability if:

- The device has been misused, abused, tampered with, or used improperly during this period.
- Malfunction results from inadequate cleaning or failure to follow the instructions.
- The equipment is operated or maintained outside the parameters indicated in the Drive DeVilbiss operating and service instructions.
- Unqualified service personnel conduct routine maintenance or servicing.
- Unauthorized parts or components (i.e., regenerated sieve material) are used to repair or alter the equipment.
- Unapproved filters are used with the unit.

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights that vary from state to state.

Ordering and Returning Parts

Drive DeVilbiss Customer Service Contact Information

Customer Service (USA): 877-224-0946

International Department: kontakt@devilbisshc.com

Ordering Non-Warranty Replacement Parts

Order non-warranty parts and literature from your Drive DeVilbiss provider. To expedite the process, be prepared to provide the following information:

- Account and ship-to numbers
- Ship-to address
- Part numbers and/or descriptions
- Quantity required
- Unit catalog number
- Unit serial number
- Hours of operation

Return and Disposal

This device may not be disposed of with household waste. After use of the device, please return the device to the provider for disposal. This device contains electrical and/or electronic components that must be recycled per EU Directive 2012/19/EU-Waste Electrical and Electronic Equipment (WEEE). Non-infectious used accessories (e.g. nasal cannula) can be disposed of as residential waste. The disposal of infectious accessories (e.g. nasal cannula from an infected user) must be made via an approved waste disposal company. Names and addresses can be obtained from the local municipality.

Store used sieve beds in sealed airtight container to minimize release of adsorbed substances. Waste must be handled in accordance with all federal, state, provincial, and local regulations. There are no known hazardous contaminants adsorbed into the molecular sieve therefore used sieve beds may be disposed of as residential waste. Never dispose by means of public sewers or drainage.

SPECIFICATIONS

PulmO2 10 Liter OXYGEN CONCENTRATOR	
Catalog Number	1060AW
Delivery Rate	1 to 10 LPM
Maximum Recommended Flow (@ nominal outlet pressures of zero and 7 kPa)**	10 LPM
Outlet Pressure	58.6 ± 3.5 kPa (8.5 ± 1.0 PSIG)
Oxygen Percentage	1-10 LPM = 93%-6%, +3%
Minimum volume fraction of delivered oxygen at the maximum rated flowrate	85% O2 @ 10 LPM
Flow Meter Accuracy	+/- 0.5 SLPM, at 10 SLPM; +0, -0.5 SLPM
5.1 ENVIROMENTAL	
Operating Atmospheric Pressure	1013 hPa to 795 hPa 0-2000 meters (0-6562 feet) Tested at nominal voltage.
Operating Temperature Range	5°C (41°F) to 40°C (104°F)
Operating Relative Humidity Range	15% to 95%, non-condensing
Power Consumption	230 VAC, 50Hz: 520 Watts Maximum
Sound Pressure Level at 3 and 10 LPM	<60 dBA
Sound Power Level at 3 and 10 LPM	<70 dBA
5.2 PHYSICAL	
Weight	25.2 kg (55.5 lbs)
Dimensions	70.5 cm H x 40.6 cm W x 35.6 cm D (27.75" H x 16" W x 14" D)
Storage Conditions	-25°C (-13°F) to 70°C (158°F), humidity range of 15% to 95% non-condensing
5.3 STANDARDS COMPLIANCE	
Approval Body and Safety Standard	TUV IEC 60601-1:2020 IEC 60601-1-6:2020 IEC 60601-1-11:2020 EN ISO 80601-2-69:2020
EMC Compliance To	EN60601-1-2:2020
5.4 ELECTRICAL	
Electrical Rating	100-240 VAC, 50 / 60Hz, 3.2 Amp
5.5 PRESSURE	
Maximum Limited Pressure under Normal condition under Single Fault condition	186 kPa (27 PSIG) 276 kPa (35 PSIG)
Operating System	Time Cycle / Pressure Swing
5.6 CLASSIFICATION	
Equipment Class and Type	Class II Equipment; Type BF Applied Part IP21 Ingress Protection

Units of measure for volumes, flows and leakages are expressed in STPD: 101,3 kPa at an operating temperature of 20 °C, dry

Measurement Uncertainty

Flow TSI 5210 flow meter +/- 1.25% of reading or 0.0005 SLPM

Pressure Heise Pressure meter (PSI) +/- 50 PPM of reading

Oxygen purity Servomex MiniMP 5200 Oxygen analyzer +/- 0.1% oxygen

Specifications subject to change without notice.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Oxygen Concentration vs Flow Rate

Flow L/m	%O ₂
10	87% - 96%
9	87% - 96%
8	87% - 95%
7	87% - 96%
6	87% - 96%
5	87% - 96%
4	87% - 96%
3	87% - 96%
2	87% - 96%
1	87% - 96%

Electromagnetic Compatibility Information

APPENDIX A

Electromagnetic Compatibility Information

IMMUNITY TEST LEVELS			
TEST DESCRIPTION	SPECIFICATION	NOTES	RESULTS
CISPR11 & IEC 60601-1-2:2014 + AMD1:2020 (EDITION 4.1)			
RF Emissions	CISPR 11	Group 1	Complies
RF Emissions	CISPR 11	Class B	Complies
Harmonic Current Emissions	EN 61000-3-2:2014	AC Input Class B	Complies
Voltage Fluctuations & Flicker Emissions	EN 61000-3-3:2013	AC Input	Complies
Electrostatic Discharge Immunity (ESD)	EN 61000-4-2:2008	±15kV (Air) ±8kV (Contact)	Complies
Radiated Electromagnetic Field Immunity	EN 61000-4-3:2006	80MHz to 2.7GHz @ 10V/m 2.7-6.0GHz @ 3V/m 80% AM @ 1kHz	Complies
Fast Transient/Burst Immunity	EN 61000-4-4:2004	±2kV for Power Lines ±1kV for I/O Lines	Complies
Surge Immunity	EN 61000-4-5:2006	±1kV differential ±2kV common	Complies
Conducted RF Immunity	EN 61000-4-6:2009	3Vrms 150kHz to 80MHz 6Vrms ISM Band frequencies between 150kHz-80MHz	Complies

Maximum length of AC line cord = 3.2m (9.8 ft)

WARNINGS

- Do not bring the device or accessories into a Magnetic Resonance (MR) environment as it may cause unacceptable risk to the patient or damage to the oxygen concentrator or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.
- Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the oxygen concentrator. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, disconnect the power cord and discontinue use. Contact your healthcare provider or clinical staff.
- This device is suitable for use in home and healthcare environments except for near active HF surgical equipment and the RF shielded room of a magnetic resonance imaging system, where the intensity of electromagnetic disturbances is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the oxygen concentrator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

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Drive DeVilbiss Healthcare Ltd.

Sidhil Business Park
Holmfield, Halifax
West Yorkshire HX2 9TN
ENGLAND
+44 (0) 845 0600 333
www.drivedevilbiss.co.uk

Drive DeVilbiss Australia Pty Limited

Building F, 2 Hudson Ave
Castle Hill, NSW 2154
AUSTRALIA
+61 02 9899 3144
www.drivedevilbiss.com

Drive DeVilbiss Healthcare France

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DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501 • USA

800-338-1988 • 814-443-4881 • www.drivemedical.com

DeVilbiss® and OSD® are registered trademarks of DeVilbiss Healthcare. Teflon® is a registered trademark of DuPont. Snoop® is a registered trademark of SWAGELOK

© 2024 DeVilbiss Healthcare LLC 10.23 All Rights Reserved.

LT-2355 Rev C

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Drive Medical GmbH & Co. KG